

# IHE Integration Statement

Product Name: ACUSON P Family  
Ultrasound System

Version:  
ACUSON P500™ – VB10

Date: October 2016

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	ACUSON P Family Ultrasound System		Oct 2016
	ACUSON P500	VB10	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
<sup>(a)</sup> Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000)	
		Vascular Ultrasound Procedure Report (TID 5100)	
Cardiology Evidence Documents	Acquisition Modality	Echo Evidence  Adult Echocardiography Procedure Report (TID 5200)	
Internet address for vendor’s IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	