

IHE Integration Statement

Product Name: ACUSON SC2000™
Volume Imaging Ultrasound
System

Version: VB20

Date: February 2017

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON SC2000 Volume Imaging Ultrasound System	VB20	Feb 2017
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cardiology			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Evidence Documents	Acquisition Modality	Echo Evidence	
	Evidence Creator (Templates: TID 5100 Vascular Ultrasound TID 5200 Echocardiography)	Echo Evidence	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org	