



INSTITUT
JULES BORDET
INSTITUUT

Routine DE Acquisitions in an Oncological Center

Patrick Emonts, MD

Medical Imaging

Institut Jules Bordet, Brussels, Belgium

Université Libre de Bruxelles (ULB)

Routine DE Acquisitions in an Oncological Center: Structure

- ♦ Introduction on CT in oncology
- ♦ Technical choices in routine acquisitions
- ♦ Clinical DE applications selection
 - ♦ VNC - iodine map use in current exams
 - ♦ Body CT examples in oncology
 - ♦ Perfusion imaging and quantification in routine pulmonary CT angiography
- ♦ Bonus (SE) Abdominal perfusion application
 - ♦ Hepatic perfusion with arteriography extraction in radioembolization planning

Intro - CT-scans Specificities in Oncology

- ◆ Screening & FU are the most usual indications
- ◆ > 80% CECT
- ◆ > 50% multiple regions CT (all body, thorax + H&N...)
- ◆ Many CT as part of clinical studies

➔ We need standardized techniques

Intro – Patients Specificities in Oncology

- ◆ Frequent examinations resulting in high cumulative radiation dose
- ◆ Longer life expectancy
- ◆ At high risk of pulmonary embolism
- ◆ Increased prevalence of renal failure

➔ Limit iodine and radiation doses

In House Routine Technical Choices





- ◆ Somatom Force® - **All CECT in dual energy**
- ◆ Common use of VNC & iodine map
- ◆ **No more non-enhanced CT before a CECT**
- ◆ Maximal use of IR algorithms and dose reduction
- ◆ Very low dose CT for pulmonary screening (DLP 20 mGy.cm)
- ◆ Reduced iodine quantity for pulmonary angiography (40ml, 16g)
- ◆ Abdominal CT: iodine volume/flow related to the patient's weight

Abdominal CT : iomeprol 400 mg/ml

Best portal phase

< 56 kg	75 ml	2.5 ml/s
56-65 kg	90 ml	3 ml/s
66-85 kg	105 ml	3.5 ml/s
86-95 kg	120 ml	4 ml/s
>95 kg	130 ml	4.5 ml/s


VNC-Iodine Map in Current Exams



Density ROI on VNC-Iodine Map Fusion

Various Clinical Applications

- ◆ Know the spontaneous density of adrenal lesion
- ◆ Be sure an hepatic lesion is enhancing
- ◆ Differentiate hematoma and enhancement
- ◆ Estimate enhancement



◆ ...
Calculated VNC density ~ 10 HU lower than measured density!

Density ROI and Iodine fraction

Potential functional studies

- ◆ Follow an hypervascular lesion
 - ◆ Chemotherapy ongoing
 - ◆ Mid-term RT
 - ◆ → early response evaluation
- ◆ Predict tumor response to anti angiogenic
- ◆ Differentiate infiltrated small nodes (5-15 mm)
- ◆ ...



Density ROI oncological examples

Case1 : FU Lung cancer post RT



Density ROI



Case1 : FU Lung cancer post RT



VNC 52.4 HU, CM -5.2 HU


Density ROI

Case1 : FU Lung cancer post RT



Density ROI


Case1 : FU Lung cancer post RT



- ◆ Conclusions:
 - ◆ Post-radiation fibrosis
 - ◆ Phantom node

Density ROI


Case 2 : FU Melanoma



- ◆ Pleural met surgery scar
- ◆ Complete remission


Density ROI

Case 2 : FU Melanoma




Density ROI

Case 2 : FU Melanoma




Density ROI

Case 2 : FU Melanoma : 4 and 8m before




Density ROI

Case 2 : FU Melanoma



- ◆ Liver
 - ◆ Low density
 - ◆ Non-enhancing
 - ◆ → Benign cyst



- ◆ Spleen
 - ◆ Low density
 - ◆ Weakly enhancing
 - ◆ Increasing size
 - ◆ → Melanoma met

Density ROI


Case 3 : FU breast cancer



- ◆ Bone, hepatic and pleural mets
- ◆ Herceptine - trastuzumab
- ◆ Complete remission

Density ROI



Case 3 : FU breast cancer



- ◆ Pleural lesions
- ◆ High density on mixed image
- ◆ Diagnosis?
 - ◆ Pleural plate
 - ◆ Met scare
 - ◆ Active met


Density ROI

Case 3 : FU breast cancer



Density ROI


Case 3 : FU breast cancer



- ◆ Very high density on VNC
→ diffusely calcified
- ◆ Substantial enhancement
→ vascularized


Density ROI

Case 3 : FU breast cancer 3y before



Density ROI

Case 3 : FU breast cancer





- ◆ Conclusions
 - ◆ Residual pleural metastases
 - ◆ Diffusely calcified
 - ◆ Latent / non-evolutive

Pulmonary CT Angiography

- ♦ DE acquisition - Systematic perfusion imaging
- ♦ DLP: 120 mGy.cm, Iodine CA: 16g
- ♦ Reduced procedure duration:
 - ♦ Exam <10 min
 - ♦ Acquisition <5 s
 - ♦ Reconstructions < 3 min
- ♦ Multiple informations obtained
 - ♦ Chest morphological CT scan
 - ♦ Angiography
 - ♦ Perfusion imaging
 - ♦ Volumes and perfusion quantification


Pulmonary CT Angiography

1 Anatomical mixed images





Pulmonary CT Angiography

2 Angiography (low KV tube)




Pulmonary CT Angiography

3 2D perfusion fusion images



Pulmonary CT Angiography

3D perfusion images



Numerous perfusion defects
< peripheral PE

Pulmonary CT Angiography

4. lung volumes & perfusion quantification

	Volume [cm3]	Full Volume [cm3]	Below Min [% of Voxels]	Above Max [% of Voxels]	Enhancement [HU]
Total	3742	3742	0	8	45+/-18
Left and Right Lung					
Right	2088	2088	0	8	42+/-15
Left	1654	1654	0	9	49+/-21
Parts of Right Lung					
Right - Upper	702	702	0	5	40+/-12
Right - Middle	691	691	0	7	37+/-13
Right - Lower	696	696	0	12	49+/-15
Parts of Left Lung					
Left - Upper	556	556	0	5	49+/-28
Left - Middle	551	551	0	8	46+/-13
Left - Lower	547	547	0	14	53+/-17

- ◆ Automatic lung segmentation & perfusion info
- ◆ Perf corresponds to average enhancement

Pulmonary CT Angiography

	A	B	C	D	E	F	G
1	Quantification des volumes et de la perfusion pulmonaire						
2							
3	Volume du poumon droit :	2088 ml		55,8%			
4	Volume du poumon gauche :	1654 ml		44,2%			
5	Volume pulmonaire total :	3742 ml		100%			
6							
7	Répartition perfusionnelle du poumon droit						
8	1/3 supérieur :	40 UH		31,7%	du poumon D	17,7%	du total
9	1/3 moyen :	37 UH		29,4%	du poumon D	16,4%	du total
10	1/3 inférieur :	49 UH		38,9%	du poumon D	21,7%	du total
11	Poumon D			100%	du poumon D	55,8%	du total
12							
13	Répartition perfusionnelle du poumon gauche						
14	1/3 supérieur :	49 UH		33,1%	du poumon G	14,6%	du total
15	1/3 moyen :	46 UH		31,1%	du poumon G	13,7%	du total
16	1/3 inférieur :	53 UH		35,8%	du poumon G	15,8%	du total
17	Poumon G			100%	du poumon G	44,2%	du total
18							
19	Poumons D+G					100%	du total

- Introduce in a table to calculate perf % of each 1/3 lung

Hepatic Perfusion CT – Arteriography Extraction

- ◆ Use in locoregional hepatic transarterial therapies
 - ◆ **Radioembolization**
 - ◆ Chemoembolization
 - ◆ Hepatic arterial infusion
- ◆ Radioembolization proceeds in 3 steps
 - ◆ **Perfusion CT with arteriography extraction**
 - ◆ Embolization simulation
 - ◆ Treatment

Hepatic Perfusion CT – Arteriography Extraction



- ◆ 1st step before invasive procedures
- ◆ Objectives:
 - ◆ Appreciate arterial part of tumoral perfusion
 - ◆ Predict the response to treatment
 - ◆ Obtain a CT hepatic arteriography to plan the selective arteriography

Hepatic Perfusion CT – Arteriography Extraction Clinical Case

- ◆ CRC with exclusive hepatic bilobar mets
- ◆ Progressive after 2 chemotherapy runs
- ◆ FDG-PET + lesions
- ◆ Full liver radioembolization decision



Hepatic Perfusion: Post Processing



- ◆ HPI

- ◆ Lesions 1 86%
- ◆ Lesion 2 80%
- ◆ Normal tissue 36%

- ◆ → Vascularisation


- ◆ Lesions 80% arterial
- ◆ Normal T 1/3 arterial
- ◆ Good candidate for radioembolization

Hepatic Perfusion: Post Processing Arteriography Extraction

- ◆ Best arterial time point selected



Hepatic Perfusion: Post Processing Arteriography Reconstructions



Conclusions

- ◆ Technical standardization prime in oncology
- ◆ Iodine and radiation doses have to be restricted for these very exposed patients
- ◆ We recommend a systematic use of dual energy acquisitions for CECT
- ◆ We recommend to execute all the body scans on the same scanner
- ◆ VNC/Iodine Map use is often helpful

Conclusions

- ◆ A pulmonary dual energy CT angiography offers simultaneously
 - ◆ An anatomical chest imaging
 - ◆ An angiography
 - ◆ A perfusion lung imaging
 - ◆ Perfusion & volumes quantification

Conclusions

- ◆ Hepatic perfusion CT in locoregional transarterial therapies planning allows
 - ◆ To appreciate the arterial part of tumoral perfusion
 - ◆ To extract a high contrast arteriography



INSTITUT
JULES BORDET
INSTITUUT

Thank you for your attention