

IHE Integration Statement

Product Name: *syngo* MI

Applications / Symbia

Version: VB20A

Date: January 2017

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	syngo MI Applications / Symbia	VB20A	January 2017
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
NM Image	Acquisition Modality	none	
	Image Display	General NM MPR	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	