



IHE Integration Statement

syngo.Breast Care

VB30

[siemens-healthineers.com/syngo.Breast Care](https://siemens-healthineers.com/syngo.Breast%20Care)

SIEMENS
Healthineers

IHE Integration Statement		Date	April 2018
Vendor	Product Name	Version	
Siemens Healthineers	syngo.Breast Care with syngo.via	VB30	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Display	No options defined	
	Evidence Creator	None	
Access to Radiology Information (ARI)	Image Display	Multiple Sources	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
	Portable Media Importer	No options defined	
	Display	No options defined	
Evidence Documents (ED)	Evidence Creator	No options defined	
Consistent Time (CT)	Time Client	None	
Mammography Image (MAMMO)	Image Display	None	
Digital Breast Tomosynthesis (DBT)	Image Display	None	
Internet address for vendor’s IHE information: www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	www.siemens.com/dicom		
Links to general information on IHE			
North America: http://ihe.net/IHE_Worldwide/	Europe: http://ihe.net/IHE_Worldwide/	Asia/Pacific: http://ihe.net/IHE_Worldwide/	

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/services/features included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

Siemens reserves the right to modify the design and specifications contained herein without prior notice.

Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible.

For recycled components we use the same extensive quality assurance measures as for factory-new components.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

syngo.Breast Care Workplace
VB30A

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone +49 9131 84-0
siemens.com/healthineers