

IHE Integration Statement

**Product Name: ACUSON
Redwood™ Ultrasound System**

Version: VA10

Date: October 2019

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON Redwood	VA10	Oct 2019
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000)	
		Vascular Ultrasound Procedure Report (TID 5100)	
		Adult Echocardiography Procedure Report (TID 5200)	
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Internet address for vendor’s IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: https://www.ihe-europe.net/	In Japan: http://www.ihe-j.org	