

IHE Integration Statement **SOMATOM go. platform**

Computed Tomography

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthcare GmbH	SOMATOM go.Up / SOMATOM go.Now SOMATOM go.ALL / SOMATOM go.Top SOMATOM go.Sim / SOMATOM go.Open Pro	VA30	11/2019
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
		MPPS - Modality group case → Siemens solution: Study Split resp. Multi Study	
	Image Display	No options defined	
	Evidence Creator	none	
Radiation Exposure Monitoring (REM)	Acquisition Modality	none	
Access to Radiology Information (ARI)	Image Display	none	
Consistent Presentation of Images (CPI)	Print Composer	No options defined	
Consistent Time (CT)	Time Client	No options defined	
Portable Data for Imaging (PDI)	Portable Media Creator	No options defined	
	Image Display	No options defined	
	Print Composer	No options defined	
Internet address for vendor's IHE information: https://www.healthcare.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation:			
HL7	https://www.healthcare.siemens.com/hl7		
DICOM	https://www.healthcare.siemens.com/dicom		
Links to general information on IHE:			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components. Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced. Caution: Federal law restricts this device to sale by or on the order of a physician.

EU Authorized Representative

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

Siemens Healthineers

Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Shanghai Medical
Equipment Ltd.
278 Zhou Zhu Road
201318 Shanghai
P. R. China