

## ACUSON Juniper ultrasound system

# Everything you need, nothing you don't

[siemens-healthineers.com/juniper](https://siemens-healthineers.com/juniper)



### Functional Performance

- Shear wave elastography imaging
- Comprehensive measurements package
- 16 transducers to cover all clinical applications
- New front-end, beamformer and back-end engines for improved contrast resolution, plunkability, sensitivity and specificity

#### 13.3" touch display

Smart UI for quicker scan settings and protocols

#### Intuitive control panel

Customizable height, angle and keys

#### Quiet, stress-free environment

Bedroom-like sound level, 40% quieter than average

#### Every detail designed to maximize productivity

Multiple storage bins available as well as space for easy installation of peripherals



#### 21.5" LED monitor

Highly adjustable for maximum operator adaptability

#### Lightweight maneuverability

One of the easiest systems to maneuver, weighing only 76 kilos (168 pounds), and with both rear and front handles

#### Purposeful cable management hooks

Several hooks available for a clean, cable free environment

#### Multiple transducer ports

Efficient workflow with 5 active transducer ports, 1 pencil cardiac port and 2 physio ports