

Compatibility table for CT

SOMATOM scanner	Software Version	syngo Virtual Cockpit Lite	syngo Virtual Cockpit
Emotion, Scope, Perspective	VC30	8.2	Not available
	VC40	8.2	Not available
	VC50	10.0	Available
Force, Flash, Drive, Edge Plus, Edge, AS	VA44	8.1	Not available
	VA48	8.2	Not available
	VA50	8.2	Not available
	VA62	8.2	Not available
	VB10	10.0	Not available
	VB20	10.0	Available
go.Now, go.Up, go.All, go.Top, X.Cite	VA10	Not available	Not available
	VA20	Not available	Not available
	VA30	10.0	Available
How to order		--> order via trial license	--> order via sVC pricebook