

IHE Integration Statement

MULTIX Impact

Strengthen your Image

siemens-healthineers.com/multix-impact



IHE Integration Statement		Date	May 2020
Vendor	Product Name	Version	
Siemens Healthineers	MULTIX Impact	VA11	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Patient Based Worklist Query	
		PPS Exception Management	
		Assisted Acquisition Protocol Setting	
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined	
Radiation Exposure Monitoring (REM)	Acquisition Modality	No options defined	
Portable Data for Imaging (PDI)	Portable Media Creator	No options defined	
Internet address for vendor’s IHE information: www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	www.siemens.com/dicom		
Links to general information on IHE			
North America: http://ihe.net/IHE_Worldwide/	Europe: http://ihe.net/IHE_Worldwide/	Asia/Pacific: http://ihe.net/IHE_Worldwide/	

MULTIX Impact is not commercially available in all countries.

Due to regulatory reasons its future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/services/features included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. For USA, federal law stipulates that this system may only be sold to a physician or by order of a physician.

The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice.

Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible.

For recycled components we use the same extensive quality assurance measures as for factory-new components.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Shanghai
Medical Equipment Ltd.
No. 278, Zhouzhu Road
Shanghai, 201318
P. R. China