

IHE Integration Statement

**Product Name: ACUSON SC2000™
Volume Imaging Ultrasound
System**

Version : VC10

Date : February 2020

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON SC2000 Volume Imaging Ultrasound System	VC10	Feb 2020
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Portable Data for Imaging	Portable Media Creator	DVD and USB	
	Portable Media Importer	DVD and USB	
	Image Display	DVD and USB	
	Report Reader	DVD and USB	
	Print Composer	DVD and USB	
Cardiology			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Evidence Documents	Acquisition Modality	Echo Evidence	
	Evidence Creator (Templates: TID 5100 Vascular Ultrasound TID 5200 Echocardiography)	Echo Evidence	
Internet address for vendor's IHE information: https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	https://www.siemens-healthineers.com/services/it-standards/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: https://www.ihe-europe.net/	In Japan: http://www.ihe-j.org	