

IHE Integration Statement


ARTIS VE30A

Product Name

ARTIS icono floor

ARTIS icono biplane

ARTIS pheno



IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthineers	ARTIS icono floor ARTIS icono biplane ARTIS pheno	VE30A	02/2023
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
	Image Display	No options defined	
	Evidence Creator	None	
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
	Image Display	No options defined	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Access to Radiology Information	Image Display	Multiple Sources	
Consistent Presentation of Images	Print Composer	None	
Consistent Time	Time Client	SNTP	
Portable Data for Imaging	Portable Media Creator	None	
	Image Display	No options defined	
	Print Composer	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Internet address for vendor's IHE information: https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise			
Links to Standards Conformance Statements for the Implementation			
HL7	n. a.		
DICOM	https://www.siemens-healthineers.com/services/it-standards/dicom		
Link to general information on IHE			
IHE	https://www.ihe.net/IHE_Worldwide		

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Manufacturer
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
