


ACUSON Sequoia

Change the way you ultrasound

siemens-healthineers.com/acuson-sequoia



Address the challenges and pressures in your practice

As a clinician, you are aware of the demands of today's healthcare environment. Variability in imaging and manual tasks add to stress and time constraints you face. Streamline your workflow, reduce variability, and alleviate physical discomfort with AI Abdomen.



Expedite your exams

Automatically recognizes and labels 17 anatomical views and calculates 12 key measurements in milliseconds.




Make excellence easy

Automations standardize imaging across users and highlights when protocols aren't met to ensure protocol-driven image capture and improve exam throughput.




Save the strain and pain

Automated labeling and measuring reduce manual movements to mitigate the fatigue and injury that can result in work-related pain.



AI Abdomen reduces hand motion 47% compared to manual.*



AI Abdomen reduces hand motion 24% compared to using a protocol.*

AI Abdomen can be used without protocol scanning to standardize imaging across users for greater consistency.

* Data on file


SIEMENS
Healthineers

Confident clinical imaging to expedite exams

Expand your expertise with AI innovations designed to improve diagnostic accuracy and patient outcomes.




AI Abdomen recognizes the right transverse kidney in this image utilizing the 5C1 transducer. Notice the system also places the appropriate body marker and transducer orientation on the image.



AI Abdomen automatically recognizes, labels, and measures anatomy, placing images accurately within the utilized protocol.



17 anatomical views for a standard abdominal exam provide a comprehensive assessment of the abdomen, covering all vital anatomical structures.



Calculates 12 key measurements in milliseconds improving accuracy and standardization across users.



For more information please contact your sales representative or visit siemens-healthineers.com/sequoia

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
22010 S.E. 51st Street
Issaquah, WA 98029, USA
Phone: 1-888-826-9702
siemens-healthineers.com/ultrasound