

Nuzzle™ Veterinary Array

Get closer to your patient™

siemens-healthineers.us/accessories



Enhanced Clarity and Diagnostic Precision

Experience unmatched clarity and detail with the Nuzzle™ veterinary array for small and large animal imaging. Our super flex veterinary MR coil achieves impressive temporal and sub-millimeter spatial resolution that sets a new clinical excellence standard.



- 80% lighter compared to conventional coils, easier to position
- One multi-purpose coil that eliminates the need to change coils
- Flexible coil provides faster throughput that may reduce length of sedation
- Superior images decrease the need for additional diagnostic procedures
- One versatile veterinary array can image soft tissue diseases, ligament and meniscal injuries, cranial and vertebral anomalies, assess bronchial plexus and muscular diseases

Strategic
Partner




X-ray Products



Improved Clinical Images

- Adaptable super-flexible 24-channel two-pad body contouring RF coil for superb high-resolution imaging of organs and neurological anomalies
- Sub-millimeter spatial resolution for demanding musculoskeletal imaging, providing detailed visualization of joints with fine tissue interfaces
- High coil density and specialized geometry maximize parallel imaging, significantly improving imaging speed



Technical Specifications

Field Strength	3.0 Tesla	1.5 Tesla
Coil	Two 12 channel pads	Two 12 channel pads
Weight	0.8 lbs. — weight with cable 3.3 lbs.	0.8 lbs. — weight with cable 3.0 lbs.
Size	15.0 in x 11.4 in.	15.0 in x 11.4 in.
Thickness	0.6 in to 1.3 in.	0.6 in to 1.3 in.
Active Area	10.5 in. x 11.4 in.	10.5 in. x 11.4 in.
Number of Elements	24-channels (12 in each pad)	24-channels (12 in each pad)
Cable Length	55 in.	53 in.
Model Number	SBA24SH30D	SBA24SH15D
Print design MRI system compatibility	Whimsical	Whimsical

[➤ Learn more](#)

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
siemens-healthineers.com

USA

Siemens Medical Solutions USA, Inc.
Healthcare
40 Liberty Boulevard
Malvern, PA 19355-9998, USA
siemens-healthineers.us

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

For USA only.