



IHE Integration Statement

syngo[®] **MR**

siemens-healthineers.com

SIEMENS
Healthineers 

IHE Integration Statement			Date	April 2022
Vendor	Product Name		Version	
Siemens Healthineers	MAGNETOM Altea	MAGNETOM Avanto Fit	syngo® MR	
	MAGNETOM Amira	MAGNETOM Semptra		
	MAGNETOM Lumina	MAGNETOM Vida Fit		
	MAGNETOM Sola	MAGNETOM Vida		
	MAGNETOM Sola Fit	MAGNETOM Free.Max		
	MAGNETOM Aera	MAGNETOM Skyra		
	MAGNETOM Prisma	MAGNETOM Skyra Fit		
	MAGNETOM Free.Star			
	This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented		Options Implemented	
Scheduled Workflow(.b)	Acquisition Modality		Broad Worklist Query	
			Patient Based Worklist Query	
			Billing and Material Management	
	Evidence Creator	None		
Image Display	No options defined			
Patient Information Reconciliation	Acquisition Modality		No options defined	
Access to Radiology Information	Image Display		None	
Portable Data for Imaging	Portable Media Creator		None	
	Portable Media Importer		No options defined	
Simple Image and Numeric Report	Report Creator		Comprehensive SR	
Invoke Image Display	Image Display		No options defined	
Consistent Time	Time Client		No options defined	
Consistent Presentation of Images	Acquisition Modality		No options defined	
	Evidence Creator		No options defined	
	Image Display		No options defined	
	Print Composer		None	
Internet address for vendor's IHE information: www.siemens-healthineers.com/services/it-standards/				
Links to Standards Conformance Statements for the Implementation				

HL7	www.siemens-healthineers.com/services/it-standards/hl7	
DICOM	www.siemens-healthineers.com/services/it-standards/dicom	
Links to general information on IHE		
North America: http://ihe.net/IHE_Worldwide/	Europe: http://ihe.net/IHE_Worldwide/	Asia/Pacific: http://ihe.net/IHE_Worldwide/

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone +49 9131 84-0
[siemens.com/healthineers](https://www.siemens.com/healthineers)