

DirectSetup Notes

Enhancing radiation therapy workflows through optimizing patient setup reporting at the CT scanner

DirectSetup Notes¹ streamlines radiotherapy patient simulation documentation by facilitating the capture of setup images and the recording of positioning accessories via tablets integrated with the CT simulator.



Simplify your workflow

This advanced tool supports radiotherapy staff to seamlessly input comments directly into customizable templates, which are designed to be versatile, adapting to a variety of treatment regions and clinical scenarios. DirectSetup Notes collects the recorded information in a comprehensive document, subsequently sending it to the Oncology Information System (OIS) for future reference and analysis.

Its user-friendly interface assists in:

- Improving efficiency in capturing patient setup images and documenting patient positioning accessories, leveraging the convenience of tablets linked to the CT simulator and its capacity to directly capture comments during patient simulation.
- Reducing numerous interactions between the CT scanner room and the acquisition console.
- The re-simulation workflow allows to easily reproduce the initial setup documentation².

The benefits of DirectSetup Notes

1. Simple

Setup information is sent as a Word document or a PDF report and can be imported² to the correct fields in OIS

2. Customized

Configure the report layout and the accessories of your own clinical practice

3. Efficient

Record your notes and patient setup images directly at the CT scanner

Enhancing your workflow in just 4 steps



- Position the patient using immobilization and additional devices as usual
- Start the simulation process



- Document setup observations directly onto the integrated tablet of the CT scanner



- Capture setup images directly using the tablet's camera



- Send the report from the CT scanner to the OIS simply and seamlessly

Siemens Healthineers Headquarters
Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Disclaimers:

¹ DirectSetup Notes is intended for digitalizing patient setup documentation to ease the RT workflow. The operation and the configuration of the application must be conducted by trained healthcare professionals.

² Requires syngo CT VB20, ARIA CORE 18.1 MR1

The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed.