

# IHE Integration Statement

**Product Name: ACUSON  
Redwood™ Ultrasound System**

**Version: VA20**

**Date: January 2021**

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON Redwood™	VA20	Jan 2021
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000)	
		Vascular Ultrasound Procedure Report (TID 5100)	
		Adult Echocardiography Procedure Report (TID 5200)	
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Internet address for vendor's IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	