

## Point of Care Informatics

# UniPOC™ System Device Interface List

Total number of interfaces: 72

Last updated: June 17, 2022

Compatible with UniPOC System v1.6 and above.

[siemens-healthineers.com/informatics](https://siemens-healthineers.com/informatics)

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Abaxis	Piccolo Xpress	Clinical chemistry	→			
Abbott Diabetes Care	FreeStyle Precision Pro (FSPP) / Wireless mode	Glucose, Ketone	↔	•	•	•
	PCx Precision	Glucose	→			
	Precision Xceed Pro (PXP)	Glucose, Ketone	↔	•	•	•
Abbott Point of Care	Afinion AS100 Analyzer	HbA1c	↔	•		•
	Afinion 2	HbA1c	↔	•		•
	Alere I	Multitest	→			
	ID Now	Multitest	↔	•	•	
	CELL-DYN EMERALD	Hematology	→			
	i-STAT 1 (via i-STAT-DE-DMI software)	Multitest	↔	•	•	
Analyticon	Urilyzer 100 / Urilyzer 100 Pro	Urinalysis	↔	•		
	Urilyzer 100 Pro V2	Urinalysis	↔	•	•	•
BD (Becton Dickinson)	BD Veritor	Multitest	→			
Boule Medical	Medonic M-series—includes M32B, M32M, M32C, M32S versions	Hematology	→			
	Medonic M-series—includes M16/M20M, M16M/M20M, M16C/M20C, M16S/M20S versions	Hematology	→			
	Swelab Alpha—includes Basic, Standard, Cap, Sampler versions	Hematology	→			
	Swelab Alpha—includes Plus Basic, Plus Standard, Plus Cap, Plus Sampler versions	Hematology	→			
Cepheid	GeneXpert PCR	Multitest	→			
	GeneXpert Omni	Multitest	↔		•	
	GeneXpert Xpress	Multitest	↔	•		•
Chirus	Dxpress Reader	Multitest	→			
CoaguSense	Coag-Sense	Coagulation	↔	•	•	
Drucker Diagnostics	QBC Star	Hematology	→			
Haemonetics	TEG 6s	Hemostasis	↔		•	

→ Unidirectional connectivity

↔ Bidirectional connectivity

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
HemoCue	HemoCue 201 DM Glucose	Glucose	↔	•	•	•
	HemoCue 201 DM Glucose RT	Glucose	↔	•	•	•
	HemoCue 201 DM Hemoglobin	Hemoglobin	↔	•	•	•
HemoSonics	Quantra <sup>NEW</sup>	Hemostasis	↔	•	•	
LumiraDx	LumiraDx	Multitest	↔	•	•	•
Medtronic	ACT Plus	Coagulation	→			
	HMS Plus	Coagulation	→			
MeMed	Key <sup>NEW</sup>	Immunoassay	→			
Nova Biomedical	NovaNet (software)	Supplier software	↔	•		
	StatSensor Creatinine (via NovaNet) StatSensor39	Creatinine	↔			
	StatStrip Glucose Meter (via NovaNet) StatStrip	Glucose	↔			
	StatStrip Lactate (via NovaNet) StatStrip32	Lactate	↔			
	StatStrip Ketone (via NovaNet) StatStrip79	Ketone	↔			
	Prime Plus (via NovaNet)	Blood Gas	↔			
	StatStrip Hemoglobin & Hematocrit (via NovaNet) StatStripHH	Hemoglobin	↔			
	StatStrip Hemoglobin & Hematocrit, Lactate (via NovaNet) StatStripHHL	Hemoglobin, Hematocrit, Lactate	↔			
Philips	BiliChek	Bilirubin	→			
Quidel	Sofia, Sofia 2	Immunoassay	→			
	Triage Meter Pro	Cardiac markers	↔	•		•
Radiometer	ABL80 Flex	Blood gas	↔		•	
	ABL80 Flex CO-OX	Blood gas	↔		•	
	ABL80 Flex CO-OX OSM	Blood gas	↔		•	
	ABL90 FLEX	Blood gas	↔		•	
	ABL90 FLEX Plus	Blood gas	↔		•	
Roche Diagnostics	ACCU-CHEK Inform II	Glucose	↔	•	•	•
	CoaguChek XS Plus	Coagulation	↔	•	•	•
	CoaguChek XS Pro	Coagulation	↔	•	•	•
	CoaguChek Pro II	Coagulation	↔	•	•	•
	cobas h232 (old)	Cardiac markers	↔	•	•	•
	cobas h232 (new)	Cardiac markers	↔	•	•	•
	cobas Liat	Multitest	↔	•		
	Urisys 1100 (from FW-version 5.x)	Urinalysis	↔	•		

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Siemens Healthcare Diagnostics	CLINITEK Advantus® Urine Chemistry Analyzer	Urinalysis	→			
	CLINITEK Status®+ Urine Chemistry Analyzer	Urinalysis	↔	•		
	DCA Vantage® Analyzer	HbA1c	↔	•		
	epoc® Blood Analysis System (Reader and Host) <b>UPDATE</b>	Multitest	↔	•	•	•
	epoc® Blood Analysis System with epoc® NXS Host (Reader and Host) <b>UPDATE</b>	Multitest	↔	•	•	•
	RAPIDPoint® 500 Blood Gas System	Blood Gas	↔	•	•	
	RAPIDPoint® 500e Blood Gas System (includes v 5.1)	Blood Gas	↔	•	•	
Sight	OLO	Hematology	→			
Sysmex	pocH-100i	Hematology	→			
Veda.Lab	EasyReaderPlus	Multitest	↔	•		
Werfen	Avoximeter 1000E	Blood Gas	→			
	Avoximeter 4000	Blood Gas	→			
	GEM Premier 4000	Blood Gas	↔	•	•	
	GEM Premier 5000	Blood Gas	↔	•	•	
	HEMOCHRON Signature Elite	Coagulation	↔	•		•
	ROTEM Delta	Hemostasis	↔		•	

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. By constantly bringing breakthrough innovations to market, we enable healthcare professionals to deliver high-quality care, leading to the best possible outcome for patients.

Our portfolio, spanning from in-vitro and in-vivo diagnostics to image-guided therapy and innovative cancer care, is crucial for clinical decision-making and treatment pathways. With our strengths in patient twinning, precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the biggest challenges in healthcare. We will continue to build on these strengths to help fight the world's most threatening diseases, improving the quality of outcomes, and enabling access to care.

We are a team of 66,000 highly dedicated employees across more than 70 countries passionately pushing the boundaries of what's possible in healthcare to help improve people's lives around the world.

CLINITEK Advantus, CLINITEK Status, DCA Vantage, epoc, RAPIDPoint, UniPOC, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local representative for availability.

---

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

**Published by**

Siemens Healthcare Diagnostics Inc.  
Point of Care Diagnostics  
2 Edgewater Drive  
Norwood, MA 02062-4637  
USA  
Phone: +1 781-269-3000