

# IHE Integration Statement

**Product Name:**  
**ACUSON Origin™ ICE 4.0**  
**Ultrasound System**

**Release: VD10**  
**Date: Jan 2026**

# IHE Integration Statement

| Vendor  | Product Name  | Version                      | Date     |
|---|---|------------------------------|----------|
| Siemens Healthineers  | ACUSON Origin™ ICE 4.0  | VD10                         | Jan 2026 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below.  |   |                              |          |
| Integration Profiles Implemented  | Actors Implemented  | Options Implemented          |          |
| <b>Radiology</b>  |   |                              |          |
| Scheduled Workflow  | Acquisition Modality  | Broad Worklist Query         |          |
|   |   | Patient Based Worklist Query |          |
| Patient Information Reconciliation  | Acquisition Modality  | None                         |          |
| Evidence Documents  | Acquisition Modality  | None                         |          |
|   | Evidence Creator<br>(Template:<br>TID 5000 OB-GYN Ultrasound)   | None                         |          |
| Portable Data for Imaging   | Portable Media Creator  | DVD and USB                  |          |
|   | Portable Media Importer   | DVD and USB                  |          |
|   | Image Display   | DVD and USB                  |          |
|   | Print Composer  | DVD and USB                  |          |
| <b>Cardiology</b>   |   |                              |          |
| Echocardiography Workflow   | Acquisition Modality  | Broad Worklist Query         |          |
|   |   | Patient Based Worklist Query |          |
| Evidence Documents  | Acquisition Modality  | Echo Evidence                |          |
|   | Evidence Creator<br>(Templates:<br>TID 5100 Vascular Ultrasound<br>TID 5200 Echocardiography<br>TID 5220: Pediatric, Fetal and<br>Congenital Cardiac Ultrasound<br>Reports) | Echo Evidence                |          |
| Portable Data for Imaging   | Portable Media Creator  | DVD and USB                  |          |
|   | Portable Media Importer   | DVD and USB                  |          |
|   | Image Display   | DVD and USB                  |          |
|   | Print Composer  | DVD and USB                  |          |
| <b>Internet address for vendor's IHE information:</b> <a href="https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise">https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise</a> |   |                              |          |
| <b>Links to Standards Conformance Statements for the Implementation</b>   |   |                              |          |
| <b>HL7</b>  | <a href="https://www.siemens-healthineers.com/services/it-standards/hl7">https://www.siemens-healthineers.com/services/it-standards/hl7</a>                                 |                              |          |
| <b>DICOM</b>  | <a href="https://www.siemens-healthineers.com/services/it-standards/dicom">https://www.siemens-healthineers.com/services/it-standards/dicom</a>                             |                              |          |

**Links to general information on IHE**

|   |  |   |
|---|--|---|
| In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>       | In Europe: <a href="https://www.ihe-europe.net/">https://www.ihe-europe.net/</a> | In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a> |
| In China: <a href="https://www.ihe-china.com">https://www.ihe-china.com</a> |  |   |