

IHE Integration Statement

Product Name: *syngo.via*

Version: VB10B

Date: 2016-04-19

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare GmbH	<i>syngo.via</i>	VB10B	April 2016
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Evidence Creator	None	
	Image Display	None	
Access to Radiology Information	Image Display	Multiple Sources	
Portable Data for Imaging	Portable Media Creator	None	
	Portable Media Importer	No options defined	
	Display	No options defined	
Evidence Documents	Evidence Creator	No options defined	
Invoke Image Display	Image Display	No options defined	
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	http://www.siemens.com/hl7		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net/IHE_Worldwide/	In Europe: http://www.ihe.net/IHE_Worldwide/	In Japan: http://www.ihe.net/IHE_Worldwide/	