

# IHE Integration Statement

**Product Name: ACUSON Juniper™  
Ultrasound System**

**Version: VB30**

**Date: January 2023**

# IHE Integration Statement

| Vendor   | Product Name   | Version   | Date         |
|--|--|---|--------------|
| Siemens Healthineers   | ACUSON Juniper™  | VB30  | January 2023 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |  |   |              |
| Integration Profiles Implemented   | Actors Implemented   | Options Implemented   |              |
| Patient Information Reconciliation   | Acquisition Modality   | None  |              |
| Scheduled Workflow   | Acquisition Modality   | Broad Worklist Query  |              |
|  |  | Patient-based Worklist Query                                      |              |
| Evidence Documents   | Acquisition Modality   | OB-GYN Ultrasound Procedure Report (TID 5000)                     |              |
|  |  | Vascular Ultrasound Procedure Report (TID 5100)                   |              |
|  |  | Adult Echocardiography Procedure Report (TID 5200)                |              |
| Echocardiography Workflow  | Acquisition Modality   | Broad Worklist Query  |              |
|  |  | Patient-based Worklist Query                                      |              |
| Internet address for vendor’s IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>   |  |   |              |
| Links to Standards Conformance Statements for the Implementation   |  |   |              |
| HL7  | n.a.   |   |              |
| DICOM  | <a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>        |   |              |
| Links to general information on IHE  |  |   |              |
| In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>  | In Europe: <a href="https://www.ihe-europe.net">https://www.ihe-europe.net</a> | In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a> |              |