

ACUSON Sequoia

Unexpected visual clarity and confidence for breast care

siemens-healthineers.co.uk/sequoia



Detect cancer like never before

We understand that early diagnosis and treatment of breast cancer require the highest quality of care, particularly when it comes to breast visualisation and penetration, especially in the case of dense breast tissue where greater clarity provides greater confidence.



Expect the unexpected

Experience resolution and penetration with our new HLX high frequency linear transducer, with exceptional resolution, deep penetration, and the ability to select UltraArt image processing filters that quickly deliver your preferred image.



Detect cancer like never before

See the difference between a stiff lesion and breast tissue to mitigate false negatives with our proprietary, next generation 2D SWE algorithm.



Minimise the manual

Needle enhancement in real time with one touch, provides crisp image details with fewer manual adjustments—InFocus imaging technology automatically focuses near and far field simultaneously.

Since 2008, worldwide breast cancer incidence has increased by more than

> 20%



Globally, breast cancer now represents

1 in 4


of all cancers in women




Everyone deserves every chance with early detection of breast cancer.

SIEMENS
Healthineers


New advances in breast visualisation deliver image quality confidence




2D Shear wave elastography allows users to easily evaluate tissue stiffness in breast masses like this fibroadenoma.




High resolution imaging can be achieved with the 15L4 linear transducer that utilises InFocus imaging technology to continuously optimise from the near to the far field.



Panoramic imaging allows the operator to record and measure structures that are larger than the transducers field of view as demonstrated in this image of a breast scar.



Quickly and accurately measure the distance, area, and circumference of a mass with AutoCalcs as demonstrated in this image of an irregular mass in the breast.



i For more information
please contact your sales
representative or visit
siemens-healthineers.
co.uk/sequoia

**Siemens Healthineers UK
Headquarters**
Siemens Healthcare Limited
Park View, Watchmoor Park
Camberley
Surrey, GU15 3YL
United Kingdom
siemens-healthineers.co.uk

Manufacturer
Siemens Medical Solutions USA, Inc.
Ultrasound
22010 S.E. 51st Street
Issaquah, WA 98029, USA
Phone: 1-888-826-9702
siemens-healthineers.com/ultrasound