



IHE Integration Statement

syngo Workflow SLR

VB36A



October 2025

SIEMENS
Healthineers

IHE Integration Statement		Date
Vendor Siemens Healthineers	Product Name <i>syngo</i> Workflow SLR	Version VB36A
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Information Reconciliation	Order Filler/Department System Scheduler	None
	Performed Procedure Step Manager	None
Scheduled Workflow	Order Filler/Department System Scheduler	None
	Performed Procedure Step Manager	None
Cross-enterprise Document Sharing for Imaging	Document Source	None
	Document Consumer	None
Cross-enterprise Document Sharing	Document Source	None
	Document Consumer	None
Internet address for vendor's IHE information: https://www.siemens-healthineers.com/en-in/services/it-standards/ihe-integrating-the-healthcare-enterprise		
Links to Standards Conformance Statements for the Implementation		
HL7	There is no HL7 conformance statement. The HL7 requirements are described in the <i>syngo</i> Workflow SLR Generic HIS Interface Mapping Guide.	
DICOM	www.siemens.com/dicom	
Links to general information on IHE: http://ihe.net/IHE_Worldwide/		

All product designations and company names are trademarks or registered trademarks of the corresponding companies.

Some of the specifications described herein may not be currently available in all countries. Please contact your local Siemens Healthineers Sales representative for the most current information.

syngo Workflow (Version SLR) is available in the U.S., Australia, New Zealand, and Canada only. Due to regional workflow requirements, the available features may vary slightly, depending on where the product is being offered.

Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice.

The original language of this document is English.

Siemens Healthineers Headquarters
Siemens Healthineers AG

Siemensstr. 3

91301 Forchheim

Germany

Phone: +49 9191 18-0

[siemens-healthineers.com](https://www.siemens-healthineers.com)

Published by ©Siemens Healthineers AG, 2025

October 2025