

IHE Integration Statement

Product Name:
ACUSON Sequoia™ Select
Ultrasound System

Version : VB10
Date : August 2023

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON Sequoia™ Select	VB10	Aug 2023
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Evidence Documents	Acquisition Modality	None	
	Evidence Creator (Template: TID 5000 OB-GYN Ultrasound)	None	
Portable Data for Imaging	Portable Media Creator	DVD and USB	
	Portable Media Importer	DVD and USB	
	Image Display	DVD and USB	
	Report Reader	DVD and USB	
	Print Composer	DVD and USB	
Cardiology			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Evidence Documents	Acquisition Modality	Echo Evidence	
	Evidence Creator (Templates: TID 5100 Vascular Ultrasound TID 5200 Echocardiography TID 5220: Pediatric, Fetal and Congenital Cardiac Ultrasound Reports)	Echo Evidence	
Internet address for vendor’s IHE information: https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise			
Links to Standards Conformance Statements for the Implementation			
HL7	N.A.		
DICOM	https://www.siemens-healthineers.com/services/it-standards/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: https://www.ihe-europe.net/	In Japan: http://www.ihe-j.org	