

IHE Integration Statement

Product Name: ACUSON Maple™
Ultrasound System

Version: VC10

Date: October 2025

IHE Integration Statement

| Vendor | Product Name | Version | Date |
|--|--|---------------------------------------|--------------|
| Siemens Healthineers | ACUSON Maple™ | VC10 | October 2025 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow.b | Acquisition Modality | NONE | |
| | | Broad Worklist Query | |
| | | Patient Based Worklist Query | |
| | | PPS Exception Management | |
| | | Assisted Acquisition Protocol Setting | |
| | Evidence Creator | NONE | |
| | | Creator Performed Procedure Step | |
| | | PPS Exception Management | |
| Evidence Documents | Acquisition Modality | NONE | |
| | Evidence Creator (Templates: TID 5000 OB-GYN Ultrasound Procedure Report, TID 5100 Vascular Ultrasound Report, TID 5200 Echocardiography Procedure Report, TID 5220 Pediatric, Fetal and Congenital Cardiac Ultrasound Reports) | NONE | |
| Portable Data for Imaging | Acquisition Modality | Broad Worklist Query | |
| | | Patient-based Worklist Query | |
| | Portable Media Creator | NONE | |
| | | DVD Media | |
| | | USB Media for PDI | |
| | Portable Media Importer | NONE | |
| | | DVD Media | |
| | | USB Media for PDI | |
| | Image Display | NONE | |
| | | DVD Media | |

| | | |
|--|---|---|
| | | USB Media for PDI |
| Internet address for vendor's IHE information: http://www.siemens-healthineers.com/ihe | | |
| Links to Standards Conformance Statements for the Implementation | | |
| HL7 | n.a. | |
| DICOM | http://www.siemens-healthineers.com/dicom | |
| Links to general information on IHE | | |
| In North America: http://www.ihe.net | In Europe: https://www.ihe-europe.net | In Japan: http://www.ihe-j.org |