



# MRI Acronyms

## Cross-Vendor Comparisons

 [siemens-healthineers.us/mri](https://siemens-healthineers.us/mri)



	Siemens Healthineers	GE	Philips	Canon	Hitachi
Basic Sequences					
Spin Echo	SE	SE	SE	SE	SE
Turbo Spin Echo/ Fast Spin Echo	TSE (Turbo Spin Echo)	FastSE (Fast Spin Echo)	TSE (Turbo Spin Echo)	FastSE (Fast Spin Echo)	FastSE (Fast Spin Echo)
Single-Shot TSE	HASTE	Single-Shot FSE	Single-Shot TSE	FASE	Single-Shot FSE
TSE with 90° Flip-Back Pulse	RESTORE	Fast Recovery FSE (FRFSE)	DRIVE	T2 Plus FSE	Driven Equilibrium
Hyperechoes	Hyperecho				
3D TSE with variable Flip Angle	SPACE	CUBE	VISTA	FASE3D mVox	
Reduced-FoV Imaging for 3D TSE	ZOOMit SPACE	HyperCube			
Number of Echoes in TSE	Turbo Factor	Echo Train Length (ETL)	Turbo Factor	Echo Train Length (ETL)	Echo Factor
Gradient Echo	GRE	GRE	Fast Field Echo (FFE)	Field Echo (FE)	GE
Spoiled Gradient Echo	FLASH	SPGR	T1-FFE	T1-FFE	RF Spoiled SARGE, RSSG
Coherent Gradient Echo	FISP	GRASS	FFE	SSFP	Rephased SARGE
Steady State Free Precession	PSIF	SSFP	T2-FFE		Time-Reversed SARGE, TRSG
True FISP	TrueFISP	FIESTA, COSMIC	Balanced FFE	True SSFP	Balanced SARGE, BASG

	Siemens Healthineers	GE	Philips	Canon	Hitachi
True FISP/Dual Excitation	CISS	FIESTA-C			Phase Balanced SARGE, PBSG
Double Echo Steady State	DESS				
Multi-Echo Data Image Combination	MEDIC	MERGE	M-FFE		ADAGE
Ultrafast Gradient Echo 2D with preparation pulse	TurboFLASH	Fast GRE, Fast SPGR	TFE	FastFE	RGE
Ultrafast Gradient Echo 3D with preparation pulse	MPRAGE	3D FGRE, 3D Fast SPGR, BRAVO	3D TFE	FastFE 3D	MPRAGE
Ultrafast Gradient Echo 3D with dual echo train	MP2RAGE			MP2RAGE	
Volume-Interpolated 3D GRE	VIBE	LAVA-XV	THRIVE	3D Quick	TIGRE
Susceptibility- Weighted Sequences	SWI	SWAN 2.0	SWIp	FSBB	BSI
Arterial Spin Labeling	ASL	ASL	ASL	ASL	ASL

	Siemens Healthineers	GE	Philips	Canon	Hitachi
Inversion Recovery	IR, Turbo IR (TIR)	IR, MPRIR, FastIR	IR-TSE	IR	IR
Short-Tau IR	TIRM, STIR	STIR	STIR	FastSTIR	STIR
Long-Tau IR	TIRM, Dark Fluid	FLAIR	FLAIR	FastFLAIR	FLAIR
Dual Inversion Recovery	DIR SPACE	CUBE DIR	Dual IR-TSE	Double IR	
True IR	TIR, True IR		Real IR	Real IR	Real-IR
Echo Planar Imaging, Diffusion	EPI	EPI	EPI	EPI	EPI
Number of Echoes in EPI	EPI Factor	ETL	EPI Factor	Echo Train Length (ETL)	Echo Factor
Diffusion-weighted Imaging	DWI	DWI	DWI	DWI	DWI
Apparent Diffusion Coefficient Map	ADC	ADC	ADC	ADC	ADC
Computed b-values	Computed b-values	MAGiC-DWI	Computed DWI	cDWI (computed DWI)	
Diffusion Tensor Imaging	DTI (Diffusion Tensor Imaging), MDDW (Multi-Directional Diffusion Imaging)	DTI (Diffusion Tensor Imaging)	DTI (Diffusion Tensor Imaging)	DTI (Diffusion Tensor Imaging)	DTI (Diffusion Tensor Imaging)
DTI Tractography (Fiber Tracking)	DTI Tractography	FiberTrak	FiberTrak	DTT (Diffusion Tensor Tractography)	DTI Tractography
High-resolution Diffusion Imaging	RESOLVE	MUSE PROPELLER DWI	DWI with segmented EPI	FASE DWI	RADAR DWI
Reduced-FoV Imaging for Diffusion and BOLD EPI	ZOOMit ZOOMit <sup>pro</sup>	FOCUS	Zoom Diffusion		
Turbo Gradient Spin Echo	TurboGSE, TGSE		GRASE	Hybrid EPI	

	Siemens Healthineers	GE	Philips	Canon	Hitachi
Applications					
Body Imaging					
Free breathing, motion-free 3D T1 imaging	StarVIBE (FREEZEit)		3D VANE XD		
Dynamic contrast- enhanced 3D with k-space re-ordering	TWIST-VIBE (FREEZEit)	DISCO			
Dynamic contrast- enhanced 3D with compressed sensing	Compressed Sensing GRASP-VIBE				
Body Diffusion	REVEAL	eDWI	DWIBS	Body Vision	(DWI)
Fat & Iron Evaluation	LiverLab	IDEAL-IQ	mDixon Quant		
MR Elastography	Elastography	MR-Touch	Elastography		
Women's Health, Men's Health					
High-resolution Bilateral Breast Imaging	VIEWS	VIBRANT-Flex	BLISS	RADIANCE	TIGRE
Comprehensive Prostate Imaging	SEEit				
MR Angiography					
Time of Flight	TOF	ToF, Inhance Inflow	TOF	TOF	TOF
Phase Contrast	Phase Contrast (PC)	Phase Contrast, Inhance Velocity	Phase Contrast (PC)	Phase Shift (PS)	Phase Contrast (PC)
Dynamic MRA with k-space Sharing	TWIST	TRICKS-XV	Keyhole (4D-TRAK)	DRKS	TRAQ
Non-contrast MR Angio, TSE-based	NATIVE-SPACE	Inhance Deltaflow	TRANCE	FBI, CIA	

	Siemens Healthineers	GE	Philips	Canon	Hitachi
Non-contrast MR Angio, TrueFISP-based	NATIVE-TrueFISP	Inhance Inflow IR	B-TRANCE	Time-SLIP	VASC
Non-contrast MR Angio, QISS	QISS				
Magnetization Transfer Contrast	MTC	MTC	MTC	SORS-STC	MTC
Ramped RF Pulse	TONE	Ramped RF	TONE	ISCE	SSP
Multi-slab acquisition	Multi-Slab	MOTSA	Multi-Chunk	Multi-Slab	Multi-Slab
Contrast Bolus Timing/Visualization	CARE Bolus	Smart Prep; Fluoro Triggered MRA	BolusTrak	Visual Prep	FLUTE
Cardiac Imaging					
High-resolution, free-breathing cardiac function	Compressed Sensing Cardiac Cine				
Comprehensive Cardiac Tool	BEAT (2D/3D)				
T1, T2, T2* Mapping	MyoMaps	CardioMaps	T1 Mapping, StarQuant	MOLLI (T1), T2* mapping (FFE2D mEcho)	
MSK Imaging					
Metal-Implant Imaging (2D correction)	WARP with VAT		O-MAR VAT	VAT	
Metal-Implant Imaging (3D correction)	Advanced WARP with SEMAC	MAVRIC SL	O-MAR XD		
Parametric Mapping	MapIt	CartiGram (T2)		T2 Mapping	T2* Relax Map

	Siemens Healthineers	GE	Philips	Canon	Hitachi
<b>Spectroscopy</b>					
Prostate Spectroscopy	3D CSI	PROSE	(Prostate Spectroscopy)	(Prostate Spectroscopy)	
Breast Spectroscopy	GRACE	BREASE	(Breast Spectroscopy)		Breast Spect
Weighted Elliptical Excitation for CSI	WET				
<b>Quiet Scanning</b>					
"Inaudible" 3D sequence with half radials	PETRA	SILENZ		mUTE 3D T1 (Pianissimo Zen)	
Quiet scanning with optimized gradient waveforms	QuietX	Silent Scan	ComforTone	(rounded gradient shapes)	
Quiet scanning with reduced slew rates (longer scan times)	Whisper sequences	Acoustic Reduction Technology (ART)	SofTone		SoftSound
<b>Motion Correction</b>					
Flow Compensation with Gradient Moment Nulling	GMR/Flow Comp	Flow Comp	Flow Comp; Flag	FC	GR
Motion Correction with Radial Blades	BLADE	PROPELLER 3.0	MultiVane XD	JET	RADAR
1D Navigators for Cardiac Imaging	1D PACE	Navigators	Navigators	Echo-Navigator	
2D Navigators for Abdominal Imaging	2D PACE			2D RMC	
3D Prospective Motion Correction for fMRI	3D PACE		PMC		
3D Retrospective Motion Correction for fMRI	3D ART	(BrainWave)	MC		

	Siemens Healthineers	GE	Philips	Canon	Hitachi
Motion Correction for 3D TSE with spiral navigators		PROMO			
Elastic 3D Motion Correction in Abdomen	DynaVIBE				
Soft-Tissue Motion Correction	BRACE				
Fat Suppression, Spatial Saturation					
Fat saturation—chemically	Fat Sat	Fat Sat/Chem Sat	SPIR	MSOFT	SINC, H-SINC
Fat saturation—chemically with adiabatic pulse	SPAIR	ASPIR	SPAIR	SPAIR	
Water Excitation	Water Excitation	Water Excitation	Proset	PASTA	Water Excitation
Dixon Fat-Water separation for TSE	Dixon TSE	IDEAL	mDixon XD, mDixon TSE	WFOP	FatSep
Dixon Fat-Water separation for 3D GRE	Dixon VIBE	LAVA-Flex	mDixon	WFS	FatSep
Spatial saturation	Sat Region	SAT	REST	Presat	Presat
Moving Sat Pulse	Tracking Sat	Walking Sat	Travel REST	Moving Presat	Sequential Pre Sat
Acceleration and Parallel Acquisition Techniques (PAT)					
PAT: Image-based Algorithm	mSENSE	ASSET	SENSE, dS SENSE	SPEEDER	RAPID
PAT: k-space-based Algorithm	GRAPPA	ARC			k-RAPID



	<b>Siemens Healthineers</b>	<b>GE</b>	<b>Philips</b>	<b>Canon</b>	<b>Hitachi</b>
<b>CAIPIRINHA for 3D sequences</b>	CAIPIRINHA				
<b>Integrated Calibration</b>	Integrated, Auto-Calibration	Self-Calibration (with ARC)			
<b>Separate Calibration</b>	Separate Calibration	(Calibration for ASSET)	(Calibration for SENSE)/CLEAR	(Calibration for SPEEDER)	(Calibration for RAPID)
<b>Multiple datasets calibrate each other</b>	Self-Calibration, T-PAT				
<b>Simultaneous Multi-Slice/Multi-Band</b>	SMS (Simultaneous Multi-Slice)	HyperBand	MultiBand SENSE	MultiBand SPEEDER	
<b>Compressed Sensing</b>	Compressed Sensing	HyperSense	Compressed SENSE	Compressed SPEEDER	IP-RAPID
<b>Parallel Transmission</b>					
<b>B1 Shimming</b>	TimTX TrueForm	MultiDrive	MultiTransmit	MultiPhase Transmit	(Quartet)
<b>Platform for advanced pTX applications, e.g., with accelerated multi-dimensional RF pulses</b>	TimTX TrueShape				
<b>Scanning with Continuous Table Move</b>	TimCT				
<b>Fast 2D GRE Localizer</b>	TimCT FastView				
<b>Contrast-enhanced Angio, 3D Coronal</b>	TimCT Angio				
<b>Multiple Contrasts, 2D Axial</b>	TimCT Onco				
<b>Patient Orientation Sequence</b>	Localizer, Scout	Localizer	Plan Scan	Locator	Scanogram
<b>Automated Slice Positioning</b>	AutoAlign	AIR x, Ready Brain	SmartExam	NeuroLine, SpineLine	

	Siemens Healthineers	GE	Philips	Canon	Hitachi
Sequence Parameters					
Repetition Time, Echo Time (in msec)	TR, TE	TR, TE	TR, TE	TR, TE	TR, TE
Inversion Time (in msec)	TI	TI	TI	TI	TI
Inter-Echo Spacing (TSE, EPI)	Echo Spacing	Echo Spacing	Echo Spacing	Echo Spacing	Inter-Echo Time (IET)
Averages	Averages	NEX	NSA	NAQ	NSA
Scan Measurement Time	Acquisition Time, TA	Acquisition Time	Acquisition Time	Acquisition Time	Scan Time
Field of View (FoV)	FoV [mm]	FoV [cm]	FoV [mm]	FoV	FoV
Rectangular FoV	FoV Phase, Rectangular FoV	Asymmetric FoV	Rectangular FoV	Rectangular FoV	Rectangular FoV
Shifting Slices Off Center	Off-center Shift	Off center FoV	Off-center FoV	Phase & Frequency Shift	Off-center FoV
Distance Between Slices	Distance Factor (% of slice thickness)	Gap	Gap	Gap	Slice Interval
Simultaneous Excitation	Simultaneous Excitation	POMP (Phase Offset Multiplanar)	Multi-Slice	QuadScan	Dual Slice
RF Excitation Pulse in Gradient Echo	Flip Angle	Flip Angle	Flip Angle	Flip Angle	Flip Angle
Bandwidth	Bandwidth [Hz/Px]	Receive Bandwidth [kHz]	Fat/Water Shift [pixel]	Bandwidth	Bandwidth
Variable Bandwidth	Optimized bandwidth	Variable Bandwidth	Optimized bandwidth	Matched Bandwidth	Variable Bandwidth

	<b>Siemens Healthineers</b>	<b>GE</b>	<b>Philips</b>	<b>Canon</b>	<b>Hitachi</b>
<b>Half Fourier Imaging</b>	Half Fourier, Phase Partial Fourier	Half NEX; fractional NEX	Half Scan	AFI	Half Scan
<b>Partial Echo</b>	Asymmetric Echo	Asymmetric Echo	Partial Echo	Matched Bandwidth	Half Echo
<b>Frequency Oversampling</b>	Oversampling	Anti-Aliasing	Frequency Oversampling	Frequency Wrap Suppression	Frequency Oversampling
<b>Phase Oversampling</b>	Phase Oversampling	No Phase Wrap	Fold-over Suppression	Phase wrap suppression	Anti-Wrap
<b>Segmented k-Space</b>	Lines/Segments	Views per segment	Views/Segment	Segments	Segments
<b>Time Delay/Block k-Space</b>	Time Delay	Intersegment Delay	TD	TD	
<b>Scan Synchronization with ECG</b>	ECG triggered	Cardiac Gated/ Triggering	ECG Triggered/ VCG	Cardiac Gated	ECG Triggered
<b>Delay after R-Wave</b>	Trigger Delay; TD	Trigger Delay; TD	Trigger Delay; TD	Trigger Delay; TD	Delay Time
<b>Respiratory Gating</b>	Respiratory Gated	Respiratory Comp	Trigger; PEAR	Respiratory Gated	MAR
<b>Multi-channel RF coil sensitivity normalization</b>	Prescan Normalize	PURE	CLEAR		NATURAL
<b>Central k-space filling arterial visualization</b>	Central Elliptical	Elliptic Centric	CENTRA		PEAKS

Product descriptions, comparisons, and specifications contained in this document are based on interpretation of available data at the time this material was collected and may require independent verification. Specifications were obtained from brochures, websites, and other independently published sources.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We're a leading medical technology company with over 120 years of experience and 18,500 patents globally. With about 50,000 dedicated colleagues in over 70 countries, we'll continue to innovate and shape the future of healthcare.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

---

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

**USA**

Siemens Medical Solutions USA, Inc.  
Healthcare  
40 Liberty Boulevard  
Malvern, PA 19355-9998, USA  
Phone: +1-888-826-9702  
siemens-healthineers.us