


IHE Integration Statement

YSIO X.free

VA10

siemens-healthineers.com



IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	YSIO X.pree	VA10	October 2020

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined
Consistent Presentation of Images (CPI)	Print Composer	No options defined
Radiation Exposure Monitoring (REM)	Acquisition Modality	none
Consistent Time (CT)	Time Client	No options defined

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.
DICOM	http://www.siemens.com/dicom

Links to general information on IHE

In North America: <http://www.ihe.net>

In Europe: <http://www.ihe-europe.net>

In Japan: <http://www.ihe-j.org>

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany