

# IHE Integration Statement

Product Name: ACUSON Juniper™  
Ultrasound System

Version: VC10

Date: October 2025

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON Juniper™	VC10	October 2025
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow.b	Acquisition Modality	NONE	
		Broad Worklist Query	
		Patient Based Worklist Query	
		PPS Exception Management	
		Assisted Acquisition Protocol Setting	
	Evidence Creator	NONE	
		Creator Performed Procedure Step	
		PPS Exception Management	
Evidence Documents	Acquisition Modality	NONE	
	Evidence Creator (Templates: TID 5000 OB-GYN Ultrasound Procedure Report, TID 5100 Vascular Ultrasound Report, TID 5200 Echocardiography Procedure Report, TID 5220 Pediatric, Fetal and Congenital Cardiac Ultrasound Reports)	NONE	
Portable Data for Imaging	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
	Portable Media Creator	NONE	
		DVD Media	
		USB Media for PDI	
	Portable Media Importer	NONE	
		DVD Media	
		USB Media for PDI	
	Image Display	NONE	

		DVD Media
		USB Media for PDI
Internet address for vendor's IHE information: <a href="http://www.siemens-healthineers.com/ihe">http://www.siemens-healthineers.com/ihe</a>		
Links to Standards Conformance Statements for the Implementation		
HL7	n.a.	
DICOM	<a href="http://www.siemens-healthineers.com/dicom">http://www.siemens-healthineers.com/dicom</a>	
Links to general information on IHE		
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="https://www.ihe-europe.net">https://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>