


IHE Integration Statement

# **syngo.via** **VB80A**



IHE Integration Statement		Date	October 2023
Vendor	Product Name	Version	
Siemens Healthineers	syngo.via	VB80A	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF.b)	Evidence Creator	None	
	Image Display	No options defined	
Access to Radiology Information (ARI)	Image Display	Multiple Sources	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
	Portable Media Importer	None	
	Display	No options defined	
Evidence Documents (ED)	Evidence Creator	No options defined	
Audit Trail and Node Authentication (ATNA) <sup>1</sup>	Secure Node	None	
	Secure Application	None	
Consistent Time (CT)	Time Client	None	
Internet address for vendor’s IHE information:			
<a href="https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise">https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	<a href="https://www.siemens-healthineers.com/services/it-standards/hl7">https://www.siemens-healthineers.com/services/it-standards/hl7</a>		
DICOM	<a href="https://www.siemens-healthineers.com/services/it-standards/dicom">https://www.siemens-healthineers.com/services/it-standards/dicom</a>		
Links to general information on IHE			
North America:	Europe:	Asia/Pacific:	
<a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>	<a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>	<a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>	

<sup>1</sup> will not support HL7 communications/services

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

---

Siemens Healthineers Headquarter  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
[siemens-healthineers.com](https://www.siemens-healthineers.com)

---

©Siemens Healthcare GmbH, 2023