


# Silhouette Adult Body Array

[siemens-healthineers.us/accessories](https://www.siemens-healthineers.us/accessories)



## Exceptional Resolution and Acceleration

Our Silhouette 24-channel MRI coil is FDA-cleared and designed to increase scan acceleration, provide exceptional image quality, and improve patient outcomes.



- Increased diagnostic accuracy: flexible array enables a close fit to maximize SNR and provide high spatial and temporal resolution
- High coil density maximizes parallel imaging capabilities to reduce imaging time
- Incredibly thin, pliable and foldable like a blanket to conform to all anatomy shapes
- Two pad design allows for versatility in positioning for better patient access, and IV- line/medical device management

Strategic  
Partner




X-ray Products

**SIEMENS**  
Healthineers

# Small Size. Mega Power.

- Flexible array enables a close fit to maximize the signal-to-noise ratio and provide high spatial and temporal resolution
- High coil density maximizes parallel imaging (ARC, ASSET) capabilities to improve imaging speed
- Can be used for any application including body, cardiac, orthopedic, neurologic, and musculoskeletal
- 2 pads each with 12 channels allows for smaller elements to be scanned with greater RF density



## Technical Specifications

Field Strength	3.0 Tesla	1.5 Tesla
Coil	Two 12 channel pads	Two 12 channel pads
Weight	0.8 lbs. — weight with cable 3.3 lbs.	0.8 lbs. — weight with cable 3.0 lbs.
Size	15.0 in x 11.4 in.	15.0 in x 11.4 in.
Thickness	0.6 in to 1.3 in.	0.6 in to 1.3 in.
Active Area	10.5 in. x 11.4 in.	10.5 in. x 11.4 in.
Number of Elements	24-channels (12 in each pad)	24-channels (12 in each pad)
Cable Length	55 in.	53 in.
Model Number	SBA24SH30C	SBA24SH15C
Print design MRI system compatibility	Cherry Blossom	Cherry Blossom

[↗ Learn more](#)

### Siemens Healthineers Headquarters

Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
siemens-healthineers.com

### USA

Siemens Medical Solutions USA, Inc.  
Healthcare  
40 Liberty Boulevard  
Malvern, PA 19355-9998, USA  
siemens-healthineers.us

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.