## **IHE Integraton Statement**

## MAX systems

**VF11** 

siemens-healthineers.com



## **IHE Integration Statement**

Vendor	Product Name	Version	Date
Siemens Healthineers	Multitom Rax Luminos dRF Max Luminos Agile Max LUMINOS Lotus Max Uroskop Omnia Max	VF11	October 2020
This product implements all tran Integration Profiles, Actors and C	sactions required in the IHE Techni Options listed below:	cal Framework to	support the IHE
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	<ul> <li>Patient Based</li> <li>Broad Worklist</li> <li>PPS Exception Management</li> <li>Billing and Ma Management</li> <li>Assisted Acqui Protocol Settir</li> </ul>	t Query terial isition
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defi	ned
Consistent Presentation of Images (CPI)	Print Composer	No options defi	ned
Radiation Exposure Monitoring (REM)	Acquisition Modality	none	
Portable Data for Imaging (PDI)	Portable Media Creator	No options defi	ned
Consistent Time (CT)	Time Client	No options defi	ned
Internet address for vendor's II	HE information: http://www.siem	ens.com/ihe	
Links to Standards Conformand	ce Statements for the Implement	ation	
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

## Links to general information on IHE

In North America: http://www.ihe.net
In Europe: http://www.ihe-europe.net

In Japan: http://www.ihe-j.org

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany Phone +49 9131 84-0 siemens-healthineers.com **Legal Manufacturer**Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany