At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

CLINITEK, CLINITEK Novus, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. Sysmex is a registered trademark of Sysmex Corporation. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany Phone: +49 9131 84-0 siemens-healthineers.com Published by

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637

USA

Phone: +1 781-269-3000





Improve reliability and productivity with the urinalysis technology you trust

The CLINITEK Novus® Automated Urine Chemistry Analyzer combines proven dry-pad urine chemistry technology with an easy-to-use cassette test format to ensure standardized test results and maximum productivity in busy laboratories. Offering a full spectrum of tests and a throughput of up to 240 samples per hour, the CLINITEK Novus analyzer streamlines workflow, simplifies routine urinalysis and STAT testing, and delivers consistent, high-quality results, even during peak periods.

Proven performance in urinalysis from the point of care to the central lab

CLINITEK® urine analyzers employ the same proven testing technology, regardless of the test format. That's why you can be sure of standardized test results in every setting, from point-of-care locations to the high-volume automated laboratory. The CLINITEK brand remains recognized the world over as the leader in reliable urine chemistry testing, meeting the needs of different patient populations, enhancing productivity, and providing broad clinical insights into patient health.

CLINITEK Novus PRO 12 cassette tests†	
Glucose	Protein
Bilirubin	Urobilinogen
Ketone	Nitrite
Blood	Leukocytes
рН	Specific Gravity
Albumin	Albumin-to- creatinine ratio



Unique "load-and-go" test cassettes

- Each cassette holds 450 test sets.
- Easy-to-use cassettes simplify analyzer operation, maintenance, and cleaning.
- Proven test pads ensure quality results with excellent sensitivity and specificity.
- Cassettes eliminate the need to orient traditional urinalysis test strips.
- RFID tags enable automated entry of lot number and expiration.

Broad menu

- CLINITEK Novus PRO 12 cassettes* for early detection of kidney disease.
- Albumin-to-creatinine (A:C) ratio results for microalbuminuria testing.



Tiltable color touchscreen

- Touchscreen provide intuitive navigation through a logical menu.
- Flexible menu options enable easy system customization, setup, maintenance, QC, and troubleshooting.

Convenient, onboard bar-code data entry

• Automated bar-code scanning minimizes data entry errors.



True calibration of all analytes

- Live calibration ensures confidence in results.
- Ready-to-use liquid calibration material enhances accuracy.



*Not available for sale in the U.S. Product availability varies by country. †Color and clarity are also analyzed.

200-sample capacity

- Load up to 200 samples at once.
- System maximizes productivity with true walkaway functionality.
- Analyzer accommodates a wide range of testing volumes and tube sizes.
- Liquid-level sensing speeds and simplifies sample preparation.
- Tube Imaging Library identifies tube and presence/ absence of cap to ensure pipette reliability.
- Rack handler is compatible with Atellica® UAS 800 Urine Sediment Analyzer* for a completely automated urinalysis system.