

# **IHE Integration Statement**

**Product Name: Biograph mCT, Biograph Vision**

**Version: VG8x**

**Date: June 2021**

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare GmbH	Biograph mCT, Biograph Vision	VG8x	06/2021
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
		Modality group case → Siemens solution: Study split resp. Multi Study	
	Image Display	No options defined in IHE	
	Evidence Creator	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined in IHE	
Radiation Exposure Monitoring (REM)	Acquisition Modality	No options defined in IHE	
	Dose Info Consumer (available as a standalone tool by CARE Analytics)	Except transaction “Store dose information”	
Access to Radiology Information (ARI)	Image Display	None	
Consistent Presentation of Images (CPI)	Print Composer	No options defined in IHE	
Consistent Time (CT)	Time Client	No options defined in IHE	
Internet address for vendor’s IHE information: <a href="https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise">https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	<a href="https://www.siemens-healthineers.com/services/it-standards/dicom">https://www.siemens-healthineers.com/services/it-standards/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>