

IHE Integration Statement

SOMATOM Force

SOMATOM Drive

SOMATOM Confidence

SOMATOM Edge Plus

SOMATOM Definition AS

SOMATOM Definition Edge

SOMATOM Definition Flash

Computed Tomography

| IHE Integration Statement | | | |
|--|---|--|---------|
| Vendor | Product Name | Version | Date |
| Siemens Healthcare GmbH | SOMATOM Force / SOMATOM Drive SOMATOM Confidence / SOMATOM Edge Plus SOMATOM Definition AS / Edge / Flash | VB20 | 05/2019 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow (SWF) | Acquisition Modality | Patient based worklist query | |
| | | Broad worklist query | |
| | | Modality group case → Siemens solution: Study Split resp. Multi Study | |
| | | Billing and Material Management | |
| | Image Display | No options defined | |
| | Evidence Creator | none | |
| Patient Information Reconciliation (PIR) | Acquisition Modality | No options defined | |
| Radiation Exposure Monitoring (REM) | Acquisition Modality | none | |
| | Dose Info Consumer (available as a standalone tool by CARE Analytics) | Except transaction “Store dose information” | |
| Access to Radiology Information (ARI) | Image Display | none | |
| Consistent Presentation of Images (CPI) | Print Composer | No options defined | |
| Consistent Time (CT) | Time Client | No options defined | |
| Portable Data for Imaging (PDI) | Portable Media Creator | No options defined | |
| | Image Display | No options defined | |
| | Print Composer | No options defined | |
| Internet address for vendor’s IHE information: https://www.healthcare.siemens.com/ihe | | | |
| Links to Standards Conformance Statements for the Implementation: | | | |
| HL7 | https://www.healthcare.siemens.com/hl7 | | |
| DICOM | https://www.healthcare.siemens.com/dicom | | |
| Links to general information on IHE: | | | |
| In North America: http://www.ihe.net | In Europe: http://www.ihe-europe.net | In Japan: http://www.ihe-j.org | |

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components. Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced. Caution: Federal law restricts this device to sale by or on the order of a physician.

Siemens Healthineers**Headquarters**

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany