

IHE Integration Statement

Product Name: ACUSON Juniper™
Ultrasound System

Version: VB11

Date: August 2022

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	ACUSON Juniper™	VB11	August 2022
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000)	
		Vascular Ultrasound Procedure Report (TID 5100)	
		Adult Echocardiography Procedure Report (TID 5200)	
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Internet address for vendor’s IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: https://www.ihe-europe.net	In Japan: http://www.ihe-j.org	