

IHE Integration Statement			
Vendor Siemens AG, Healthcare Sector	Product Name <i>syngo.via</i>	Version VA10A	Date Apr. 2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Display Evidence Creator	None None	
Consistent Presentation of Images (CPI)	Image Archive Image Manager	None None	
Evidence Documents (ED)	Image Archive Image Manager	None None	
Access to Radiology Information (ARI)	Image Display	Multiple Sources	
Image Fusion (FUS)	Image Archive Image Manager	None None	
Nuclear Medicine Imaging (NMI)	Image Archive Image Manager	None None	
Key Image Note (KIN)	Image Archive Image Manager	None	
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Display	None	
Consistent Time (CT)	Time Client	None	

Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	