

Point of Care Informatics

POCcelerator® System Device Interface List

Total number of interfaces: 214

Last updated: January 27, 2022

siemens-healthineers.com/informatics

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
4S Information Systems	DAWN	Supplier software	→			
77 Elektronika	DocUReader 2 Pro	Urinalysis	↔	•	•	•
Abaxis	Piccolo Xpress	Clinical chemistry	→			
Abbott Diabetes Care	FreeStyle Precision Pro (FSPP) UPDATE	Glucose, Ketone	↔	•	•	•
	FreeStyle Precision Pro (FSPP) Wireless mode UPDATE	Glucose, Ketone	↔	•	•	•
	PCx Precision	Glucose	→			
	Precision Xceed Pro (PXP)	Glucose, Ketone	↔	•	•	•
Abbott Point of Care	Afinion AS100 Analyzer	HbA1c	↔	•		•
	Afinion 2	HbA1c	↔	•		•
	Alere i	Multitest	→			
	Alere Reader	Multitest	→			
	CELL-DYN Emerald	Hematology	→			
	Cholestech LDX	Multitest	→			
	ID NOW (ASTM)	Multitest	→			
	ID NOW (POCT1-A)	Multitest	↔	•	•	•
	Pima	HIV	→			
	i-STAT 1 (via i-STAT-DE-DMI software) UPDATE	Multitest	↔	•	•	•
	i-STAT Alinity	Multitest	↔	•	•	•
Acon Laboratories	Mission U120	Urinalysis	→			
Aidian (former Orion Diagnostica)	QuikRead 101	CRP	→			
	QuikRead go (ASTM)	CRP	→			
	QuikRead go (POCT1-A2) UPDATE	CRP	↔	•	•	•
A. MENARINI Diagnostics	Zenit RA	Immunoassay	→			
Analyticon	CombiScan 100	Urinalysis	→			
	Urilyzer 100	Urinalysis	↔	•		
	Urilyzer 100 Pro	Urinalysis	↔	•	•	•
Arkray	Adams HA-8160	HbA1c	→			

→ Unidirectional connectivity

↔ Bidirectional connectivity

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Beckman Coulter	Coulter Ac-T diff	Hematology (no QC!)	→			
	Coulter Ac-T diff 2	Hematology	→			
	Coulter Ac-T 8	Hematology	→			
	Coulter Ac-T 10	Hematology	→			
	DxH 500	Hematology	→			
	DxH 520	Hematology	→			
Becton Dickinson	BD Veritor Plus (via BD Cloud Connectivity Solution)	Multitest	→			
Biomérieux	mini VIDAS	Multitest	→			
Bio-Rad	Bio-Rad Unity Real Time (URT)	Supplier software	→			
biosurfit	Spinit	Multitest	→			
Boditech	AFIAS-1	Immunoassay	→			
	AFIAS-6	Immunoassay	→			
Boule Medical	Medonic M-series M16/M20	Hematology	→			
	Medonic M-series M16C/M20C	Hematology	→			
	Medonic M-series M16M/M20M	Hematology	→			
	Medonic M-series M16S/M20S	Hematology	→			
	Medonic M-series M32B	Hematology	→			
	Medonic M-series M32C	Hematology	→			
	Medonic M-series M32M	Hematology	→			
	Medonic M-series M32S	Hematology	→			
	Swelab Alfa Basic / Plus Basic	Hematology	→			
	Swelab Alfa Standard / Plus Standard	Hematology	→			
	Swelab Alfa Cap / Plus Cap	Hematology	→			
	Swelab Alfa Sampler / Plus Sampler	Hematology	→			
Caretaker Medical	Caretaker4	Vital Signs Monitor	→			
Cepheid	GeneXpert PC and Omni ^{UPDATE}	Multitest	→			
	GeneXpert Infinity ^{UPDATE}	Multitest	→			
	GeneXpert Xpress ^{UPDATE}	Multitest	↔	•		•
Chirus	DXpress Reader	Multitest	→			
CoaguSense	Coag-Sense PT2 ^{NEW}	Coagulation	↔	•	•	
concile	Omega 100 (Ω100)	Immunoassay	→			
Diagnostics For The Real World (DRW)	Samba II ^{NEW}	Multitest	↔		•	
Diatron	Diatron Aquila	Hematology	→			
DnaNudge	NudgeBox	Immunoassay	↔	•		
Dräger	Jaundice Meter JM-105	Bilirubin	→			
Drucker Diagnostics	QBC Star	Hematology	→			
EKF Diagnostics	Quo-Lab (via CIB)	HbA1c	↔	•	•	•
	Quo-Test (via CIB)	HbA1c	↔	•	•	•
Enigma Diagnostics	Enigma MiniLab	Multitest	→			
ET Healthcare	ET Pylon	Multitest	→			

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Eurolyser Diagnostica	Eurolyser CUBE	Multitest	→			
	Eurolyser CUBE-S	Multitest	→			
	Eurolyser smart	Multitest	→			
	Eurolyser solo	Multitest	→			
Fuji	Fuji Dri-Chem 4000 UPDATE	Clinical chemistry	→			
	Fuji Dri-Chem NX10N UPDATE	Clinical chemistry	→			
	Fuji Dri-Chem NX500 UPDATE	Clinical chemistry	→			
	Fuji Dri-Chem NX600 NEW	Clinical chemistry	→			
	Fuji Dri-Chem NX700 UPDATE	Clinical chemistry	→			
GenePOC Diagnostics	Revogene	Multitest	↔	•		•
Haemonetics	ClotPro	Hemostasis	↔		•	
	TEG 6s (via TEGManager/ HaemoCommunicator)	Hemostasis	↔		•	
Helena Laboratories	Cascade Abrazo	Coagulation	↔			•
	Cascade POC	Coagulation	↔	•		
HemoCue	Glucose 201 DM	Glucose	↔	•	•	•
	Glucose 201 DM RT	Glucose	↔	•	•	•
	Hb 201 DM	Hemoglobin	↔	•	•	•
	HbA1c 501	Hemoglobin	→			
	WBC Diff	Hematology	→			
HemoSonics	Quantra	Coagulation	↔	•	•	
HLS Therapeutics	CSAN Pronto	Hematology	→			
Hologic	PeriLynx	Immunoassay	→			
	Rapid fFN	Immunoassay	→			
HORIBA	ABX Micros 60	Hematology	→			
	ABX Micros ES 60	Hematology	→			
	ABX Microsemi CRP	Hematology	→			
	Yumizen H500	Hematology	→			
Lite-On Technology Corp.	skyla HB1 (POCT1-A2)	Clinical chemistry	↔	•		
	skyla HB1 (ASTM)	Clinical chemistry	→			
	skyla Hi	HbA1c	↔	•		
LumiraDx	LumiraDx (Instrument and Host) UPDATE	Multitest	↔	•	•	•
Macherey-Nagel	Uryxson Relax (Connect)	Urinalysis	↔		•	
Medtronic	ACT Plus	Coagulation	→			
	HMS Plus	Coagulation	→			
Mitsubishi	PATHFATH	Immunoassay	→			
nal von minden (BioSciTec and Mediwiss)	Rapid Slide Scanner (RSS) III	Drug-Screen	→			
Nihon Kohden	Celltac Alpha (MEK 6500/6510)	Hematology	→			
	Celltac Chemi (CHM-4100)	Hematology	→			
Norma Diagnostika	Norma Icon 3	Hematology	→			
	Norma Icon 5	Hematology	→			

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Nova Biomedical	NovaNet*	Supplier software	↔	•	•	
	Stat Profile pHox Ultra (via NovaNet)	Blood gas	→			
	Stat Profile Prime (via NovaNet)	Blood gas	→			
	Stat Profile Prime Plus (via NovaNet)	Blood gas	↔	•	•	
	StatSensor Creatinine (via NovaNet)	Creatinine	↔	•	•	
	StatStrip Glucose (via NovaNet)	Glucose	↔	•	•	
	StatStrip Glucose/Ketone (via NovaNet)	Glucose, Ketone	↔	•	•	
	StatStrip Hb/Hct (via NovaNet)	Hemoglobin	↔	•	•	
	StatStrip HHL (via NovaNet)	Hemoglobin/ Hematocrit/Lactate	↔	•	•	
	StatStrip Lactate (via NovaNet)	Lactate	↔	•		
	Stat Profile CCX	Blood gas	↔	•	•	
	Stat Profile pHox / pHox plus	Blood gas	↔		•	
	Stat Profile pHox Ultra	Blood gas	→			
	Stat Profile Prime	Blood gas	→			
	StatSensor Creatinine	Creatinine	↔	•	•	•
	StatStrip Glucose	Glucose	↔	•	•	•
	StatStrip Glucose/Ketone	Glucose, Ketone	↔	•	•	•
	StatStrip Hb/Hct	Hemoglobin	↔	•	•	•
	StatStrip Lactate	Lactate	↔	•	•	•
OPTI Medical	OPTI CCA-TS2	Blood gas	↔	•	•	
Philips	BiliChek	Bilirubin	→			
PixCell Medical	HemoScreen	Hematology	→			
Quidel	Sofia	Immunoassay	→			
	Sofia 2	Immunoassay	→			
	Triage Meter Plus	Cardiac markers	↔	•		
	Triage Meter Pro	Cardiac markers	↔	•		
Radiometer	ABL7xx series	Blood gas	↔		•	
	ABL8xx series	Blood gas	↔		•	
	ABL80 FLEX series	Blood gas	↔		•	
	ABL90 FLEX series	Blood gas	↔		•	
	AQT90 FLEX	Multitest	↔		•	
	AQURE	Supplier software	↔		•	
	NPT7	Blood gas	→			
	RADIANCE	Supplier software	↔		•	
Response Biomedical	RAMP 200	Multitest	↔	•		•

*NovaNet is needed for the StatStrip-Family connections with a firmware version greater than 4.17.8.47.
Other devices can be connected directly or via NovaNet

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Roche Diagnostics	ACCU-CHEK Inform II UPDATE	Glucose	↔	•	•	•
	CARDIAC reader	Cardiac markers	↔			•
	CoaguChek Pro II UPDATE	Coagulation	↔	•	•	•
	CoaguChek XS Plus	Coagulation	↔	•	•	•
	CoaguChek XS Pro	Coagulation	↔	•	•	•
	cobas b 101	Hemoglobin	↔	•	•	•
	cobas b 123	Blood gas	↔	•	•	
	cobas c 111	HbA1c	→			
	cobas e 411	Immunoassay	→			
	cobas h 232 UPDATE	Cardiac markers	↔	•	•	•
	cobas 6000	Clinical chemistry	→			
	cobas IT 1000	Supplier software	↔	•	•	
	cobas Liat UPDATE	Multitest	↔	•		
	Reflotron Plus	Clinical chemistry	→			
	Reflotron Sprint	Clinical chemistry	→			
	Urisys 1100 FW5X	Urinalysis	↔	•		
Samsung	LABGEO HC10	Hematology	↔	•		
	LABGEO IB10	Immunoassay	↔	•		
	LABGEO PT10	Clinical chemistry	↔	•		
SD Biosensor	Standard F200	Multitest	→			
	Standard F2400	Multitest	→			
Siemens Healthcare Diagnostics	ADVIA 360® Hematology System	Hematology	→			
	ADVIA Centaur® CP/XP Immunoassay System	Immunoassay	→			
	Atellica® DCA Analyzer (Engine and IDM) NEW	Hemoglobin	↔	•	•	•
	Atellica® 1500 Automated Urinalysis System	Urinalysis	↔			•
	Atellica® UAS 800 Automated Urinalysis System	Urinalysis	↔			•
	Atellica® VTLi Immunoassay System UPDATE	Immunoassay	↔	•		
	CLINITEK Advantus® Urine Chemistry Analyzer	Urinalysis	→			
	CLINITEK Status® Analyzer	Urinalysis	→			
	CLINITEK Status®+ Urine Chemistry Analyzer UPDATE	Urinalysis	↔	•		•
	DCA Vantage® Analyzer UPDATE	HbA1c	↔	•		•
	epoc® Blood Analysis System (Reader and Host) UPDATE	Multitest	↔	•	•	•
	epoc® NXS Blood Analysis System (Reader and Host) UPDATE	Multitest	↔	•	•	•
	RAPIDComm® Data Management System	Supplier software	↔	•	•	
	RAPIDLab® 1200 series Blood Gas Systems	Blood gas	↔	•	•	•
	RAPIDLab® 348EX Blood Gas System	Blood gas	→			
	RAPIDPoint® 500 Blood Gas System	Blood gas	↔	•	•	•
	RAPIDPoint® 500e Blood Gas System UPDATE	Blood gas	↔	•	•	•
	Stratus® CS Acute Care™ Diagnostic System	Cardiac markers	→			
	Stratus® CS 200 Acute Care™ Diagnostic System	Cardiac markers	↔	•		
	Xprecia Stride™ (HID) Coagulation Analyzer	Coagulation	↔	•		•

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Sight Diagnostics	OLO	Hematology	→			
Sphere Medical	Proxima	Blood gas	→			
Spindiag	Rhonda Player ^{NEW}	Multitest	↔	•		
Stat-Diagnostica	DiagCORE	Multitest	→			
Sysmex	Alifax Roller	Hematology	→			
	CA-560	Coagulation	→			
	CA-600	Coagulation	→			
	poch-100i	Hematology	→			
	XN-350 (via SUIT Manager)	Hematology	→			
	XN-450 (via SUIT Manager)	Hematology	→			
	XN-550 (via SUIT Manager)	Hematology	→			
	XP-300	Hematology	→			
	XS-800i (via SUIT Manager)	Hematology	→			
Thermo Fisher Scientific	B-R-A-H-M-S Direct Reader	Multitest	→			
Veda.Lab	Easy Reader+	Multitest	↔	•		
Welch Allyn	Connex Spot Monitor	Vital Signs Monitor	→			
Werfen	Avoximeter 1000E	Blood gas	→			
	Avoximeter 4000	Blood gas	→			
	BIO-Flash ^{UPDATE}	Immunoassay	→			
	GEM Premier 3000	Blood gas	↔		•	
	GEM Premier 3500	Blood gas	↔		•	
	GEM Premier 4000	Blood gas	↔	•	•	
	GEM Premier 5000	Blood gas	↔	•	•	
	GEMweb Plus	Supplier software	↔	•	•	
	HEMOCHRON Jr. Signature	Coagulation	→			
	HEMOCHRON Response	Coagulation	→			
	HEMOCHRON Signature Elite	Coagulation	↔	•		•
	ProTime InRhythm	Coagulation	→			
	ROTEM delta (direct or via Rotem live)	Hemostasis	↔		•	
	ROTEM gamma (direct or via Rotem live)	Hemostasis	↔		•	
	ROTEM sigma (direct or via Rotem live)	Hemostasis	↔		•	

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Acute Care, ADVIA, ADVIA Centaur, Atellica, CLINITEK Advantus, CLINITEK Status, DCA Vantage, epoc, POCcelerator, RAPIDComm, RAPIDLab, RAPIDPoint, Stratus, Xprecia Stride, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000