

## IHE Integration Statement

# LUMINOS Q.namix VB10E (and higher)

### Product Name

---

LUMINOS Q.namix T

---

LUMINOS Q.namix R

---



IHE Integration Statement		Date      November 2025
Vendor Siemens Healthineers	Product Name LUMINOS Q.namix	Version VB10E (and higher)
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF.b)	Acquisition Modality	Broad Worklist Query
		Billing and Material Management
		PPS Exception Management
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined
Consistent Presentation of Images (CPI)	Print Composer	No options defined
Radiation Exposure Monitoring (REM)	Acquisition Modality	None
Portable Data for Imaging (PDI)	Portable Media Creator	No options defined
Consistent Time (CT)	Time Client	No options defined
Internet address for vendor's IHE information: <a href="http://www.siemens-healthineers.com/ihe">www.siemens-healthineers.com/ihe</a>		
Links to Standards Conformance Statements for the Implementation		
HL7	<a href="http://www.siemens-healthineers.com/hl7">www.siemens-healthineers.com/hl7</a>	
DICOM	<a href="http://www.siemens-healthineers.com/dicom">www.siemens-healthineers.com/dicom</a>	
Links to general information on IHE		
North America: <a href="http://ihe.net/IHE_Worldwide">ihe.net/IHE_Worldwide</a>	Europe: <a href="http://ihe-europe.net">ihe-europe.net</a>	Asia/Pacific: <a href="http://ihe.net/IHE_Worldwide">ihe.net/IHE_Worldwide</a>

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

---

Siemens Healthineers Headquarters  
Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
Phone: +49 9191 18-0  
siemens-healthineers.com