

## IHE Integration Statement

# ***MAMMOMAT B.brilliant VA10***



IHE Integration Statement		Date    July 2023
Vendor Siemens Healthineers	Product Name MAMMOMAT B.brilliant	Version VA10
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF / SWF.b)	Acquisition Modality	Broad Worklist Query
		Patient Based Worklist Query
		Assisted Acquisition Protocol Setting
		PPS Exception Management
	Image Display	No options defined
Patient Information Reconciliation (PIR)	Acquisition Modality	None
Mammography Acquisition Workflow (MAWF)	Acquisition Modality	No options defined
Mammography Image (MAMMO)	Acquisition Modality	None
Digital Breast Tomosynthesis (DBT)	Acquisition Modality	None
Radiation Exposure Monitoring (REM)	Acquisition Modality	No options defined
Access to Radiology Information (ARI)	Image Display	None
Portable Data for Imaging (PDI)	Portable Media Creator	DVD Media
	Image Display	DVD Media
	Print Composer	DVD Media
Consistent Presentation of Images (CPI)	Print Composer	None
Consistent Time (CT)	Time Client	None
Internet address for vendor’s IHE information: <a href="http://www.siemens-healthineers.com/ihe">www.siemens-healthineers.com/ihe</a>		
Links to Standards Conformance Statements for the Implementation		
HL7	<a href="http://www.siemens-healthineers.com/hl7">www.siemens-healthineers.com/hl7</a>	
DICOM	<a href="http://www.siemens-healthineers.com/dicom">www.siemens-healthineers.com/dicom</a>	

Links to general information on IHE		
North America: <a href="http://ihe.net/IHE_Worldwide">http://ihe.net/IHE_Worldwide</a>	Europe: <a href="http://ihe.net/IHE_Worldwide">http://ihe.net/IHE_Worldwide</a>	Asia/Pacific: <a href="http://ihe.net/IHE_Worldwide">http://ihe.net/IHE_Worldwide</a>

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

---

#### **Siemens Healthineers Headquarters**

Siemens Healthcare GmbH

SHS DI XP

Henkestr. 127

91052 Erlangen, Germany

Phone: +49 9131 84-0

siemens-healthineers.com

---