

# IHE Integration Statement

## MULTIX Impact C

VA20, VA21 and higher

[siemens-healthineers.com](https://www.siemens-healthineers.com)



# IHE Integration Statement

| Vendor   | Product Name  | Version  | Date       |
|--|---|--|------------|
| Siemens Healthineers   | MULTIX Impact C   | VA20 and<br>VA21   | March 2022 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |   |  |            |
| Integration Profiles Implemented   | Actors Implemented  | Options Implemented  |            |
| Scheduled Workflow (SWF)   | Acquisition Modality  | Patient Based Worklist Query<br>PPS Exception Management<br>Assisited Acquisition Protocol Setting |            |
| Patient Information Reconciliation (PIR)   | Acquisition Modality  | No options defined   |            |
| Radiation Exposure Monitoring (REM)  | Acquisition Modality  | No options defined   |            |
| Portable Data for Imaging (PDI)  | Acquisition Modality  | No options defined   |            |
| Internet address for vendor's IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>   |   |  |            |
| Links to Standards Conformance Statements for the Implementation   |   |  |            |
| HL7  | n.a.  |  |            |
| DICOM  | <a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a> |  |            |

## **Links to general information on IHE**

In North America: <http://www.ihe.net>

In Europe: <http://www.ihe-europe.net>

In Japan: <http://www.ihe-j.org>

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

---

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone +49 9131 84-0  
siemens-healthineers.com

**Manufacturer**

Siemens Shanghai  
Medical Equipment Ltd.  
No. 278, Zhouzhu Road  
Shanghai, 201318  
P. R. China