


IHE Integration Statement

syngo Dynamics VA50

syngo Dynamics

VA50



IHE Integration Statement		Date
		June 2026
Vendor Siemens Healthineers	Product Name syngo Dynamics	Version VA50
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Displayable Reports (DRPT – CARD)	Report Creator	No options defined
Evidence Documents (ED – CARD)	Evidence Creator	<ul style="list-style-type: none"> ▪ Echo Evidence <ul style="list-style-type: none"> ○ Vascular Ultrasound (5100) ○ Echocardiography (5200)
Cath Workflow	Image Display	No options defined
Echo Workflow	Image Display	No options defined
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	No options defined
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer	No options defined
Internet address for vendor’s IHE information: www.siemens-healthineers.com/ihe		
Links to Standards Conformance Statements for the Implementation		
HL7	www.siemens-healthineers.com/hl7	
DICOM	www.siemens-healthineers.com/dicom	
Links to general information on IHE		
North America: https://ihe.net/IHE_Worldwide/	Europe: https://ihe-europe.net/	Asia/Pacific: https://ihe.net/IHE_Worldwide/

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Legal Manufacturer

Siemens Healthineers AG
Digital & Automation
Siemensstr. 1
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Published by Siemens Healthineers AG- 062026 online only- ©Siemens Healthineers AG, 2026