


# Chek-Stix Urinalysis Liquid Quality Control Products

siemens-healthineers.com



## Quality control made effortless

As a leader in urinalysis, Siemens Healthineers strives to deliver quality solutions you can trust. Healthcare and laboratory professionals rely on Siemens Healthineers to ensure consistent, accurate results and help streamline operations. Building on that legacy of leadership and trust, Siemens Healthineers is proud to deliver the latest innovation in its urinalysis portfolio—introducing Chek-Stix® Liquid Quality Control solutions.



---

### Key benefits:



Liquid-ready solution means less time prepping quality control solution and more time running patient samples.



Barcode-enabled for quick, accurate data input.



Supports our broad urine chemistry and sediment parameter menu.



Validated for use with all Siemens Healthineers urinalysis devices and reagents.



22 months' shelf life from manufacture date.

---

# Chek-Stix Liquid Control comes in two configurations to meet each lab's unique needs

## Chek-Stix Liquid QC Kit + Sediment



### Analyte Profile

#### Chemistry

pH	Urobilinogen
Protein	Specific Gravity
Glucose	Leukocyte Esterase
Ketone	Creatinine
Bilirubin	Albumin
Blood (Hemoglobin)	A:C (Albumin-to-Creatinine Ratio)
Nitrite	P:C (Protein-to-Creatinine Ratio)

#### Sediment

RBC  
WBC

### System Compatibility

Atellica® 1500 Automated Urinalysis System\*



### Ordering Information

SMN	Each kit contains
11561696	<ul style="list-style-type: none"> <li>Level 1 (Negative) Control 1 x 120 mL</li> <li>Level 2 (Positive) Control 1 x 120 mL</li> </ul>

\*Not available for sale in the U.S.

## Chek-Stix Liquid QC Kit



### Analyte Profile

#### Chemistry

pH	Urobilinogen
Protein	Specific Gravity
Glucose	Leukocyte Esterase
Ketone	Creatinine
Bilirubin	Albumin
Blood (Hemoglobin)	A:C (Albumin-to-Creatinine Ratio)
Nitrite	P:C (Protein-to-Creatinine Ratio)

### System Compatibility

CLINITEK Novus® Automated Urine Chemistry Analyzer



### Ordering Information

SMN	Each kit contains
11561695	<ul style="list-style-type: none"> <li>Level 1 (Negative) Control 1 x 120 mL</li> <li>Level 2 (Positive) Control 1 x 120 mL</li> </ul>

Atellica, Chek-Stix, CLINITEK Novus, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

### Siemens Healthineers Headquarters

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

### Published by

Siemens Healthcare Diagnostics Inc.  
Point of Care Diagnostics  
2 Edgewater Drive  
Norwood, MA 02062-4637  
USA  
Phone: +1 781-269-3000