


Award-Winning Featherlight Snuggle Pediatric Body Array

siemens-healthineers.us/accessories



Exceptional Resolution and Acceleration

The first child-centered design for pediatric patients. Snuggle our huggable 24-channel MRI coil is FDA cleared and is designed to increase scan acceleration, provide exceptional image quality, and improve patient outcomes.



- Increased diagnostic accuracy: flexible array enables a close fit to maximize SNR and provide high spatial and temporal resolution
- High coil density maximizes parallel imaging capabilities to reduce imaging time
- Soft and ultra-light weight coil increases comfort and cooperation. Space dinosaur themed graphics and soft covers provides an engaging, friendly experience for the patient
- Separate pads improve patient access, setup, and MR suite workflow

Strategic
Partner



X-ray Products



Kid Sized. Mega Power.

- Flexible array enables a close fit to maximize the signal-to-noise ratio and provide high spatial and temporal resolution.
- High coil density maximizes parallel imaging capabilities to improve imaging speed.
- Coil elements are sized for pediatric patients and are without hard nodules so patients can lie on the array comfortably.



Patient Size Matters

[➤ Learn more](#)

Technical Specifications

Field Strength	3.0 Tesla	1.5 Tesla
Coil	Two 12 channel pads	Two 12 channel pads
Weight	0.8 lbs. — weight with cable 3.3 lbs.	0.8 lbs. — weight with cable 3.0 lbs.
Size	15.0 in x 11.4 in.	15.0 in x 11.4 in.
Thickness	0.6 in to 1.3 in.	0.6 in to 1.3 in.
Active Area	10.5 in. x 11.4 in.	10.5 in. x 11.4 in.
Number of Elements	24-channels (12 in each pad)	24-channels (12 in each pad)
Cable Length	55 in.	53 in.
Model Number	SBA24SH30B	SBA24SH15B
Print design MRI system compatibility	Space Dinos	Space Dinos

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
siemens-healthineers.com

USA

Siemens Medical Solutions USA, Inc.
Healthcare
40 Liberty Boulevard
Malvern, PA 19355-9998, USA
siemens-healthineers.us

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.