

IHE Integration Statement

Product Name: Biograph Vision Quadra

Version: VR20A

Date: February 2024

This document is same valid for future versions as long as no other document exists.

IHE Integration Statement

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthineers AG	Biograph Vision Quadra	VR20A	02/2024
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
		Modality group case → Siemens solution: Study split resp. Multi Study	
	Image Display	No options defined in IHE	
	Evidence Creator	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined in IHE	
Radiation Exposure Monitoring (REM)	Acquisition Modality	No options defined in IHE	
Access to Radiology Information (ARI)	Image Display	None	
Consistent Presentation of Images (CPI)	Print Composer	No options defined in IHE	
Consistent Time (CT)	Time Client	No options defined in IHE	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
	Image Display	None	
	Print Composer	None	
Internet address for vendor's IHE information: https://www.siemens-healthineers.com/services/it-standards/ihe-integrating-the-healthcare-enterprise			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	https://www.siemens-healthineers.com/services/it-standards/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net		In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org