

# IHE Integration Statement

## *syngo* MI Applications VB22

### Product Name

---

Symbia E, S, T

Symbia Evo, Intevo

Symbia.net

---

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthineers	Symbia E, S, T, Evo, Intevo; Symbia.net	VB22	Feb 2026
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
NM Image	Acquisition Modality	none	
	Image Display	General NM MPR	
Internet address for vendor's IHE information: <a href="http://www.siemens-healthineers.com/ihe">http://www.siemens-healthineers.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="https://www.siemens-healthineers.com/hl7">https://www.siemens-healthineers.com/hl7</a>		
<b>DICOM</b>	<a href="http://www.siemens-healthineers.com/dicom">http://www.siemens-healthineers.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	