

# IHE Integration Statement

## *syngo CT VA40A*

### Product Name

---

SOMATOM go.Platform  
SOMATOM X.cite & X.ceed  
NAEOTOM Alpha

---



<b>IHE Integration Statement</b>		<b>Date</b>	7/2021
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Siemens Healthineers	SOMATOM go.Now / go.Up SOMATOM go.All / go.Top SOMATOM go.Sim / go.Open Pro SOMATOM X.cite/X.ceed NAEOTOM Alpha	VA40	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Acquisition Modality	Patient based worklist query Broad worklist query MPPS - Modality group case → Siemens solution: Multi Study	
	Image Display	No options defined	
	Evidence Creator	none	
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined	
Radiation Exposure Monitoring (REM)	Acquisition Modality	none	
Access to Radiology Information (ARI)	Image Display	none	
Consistent Presentation of Images (CPI)	Print Composer	No options defined	
Consistent Time (CT)	Time Client	No options defined	
Portable Data for Imaging (PDI)	Portable Media Creator	No options defined	
	Image Display	No options defined	
	Print Composer	No options defined	
<b>Internet address for vendor's IHE information:</b> <a href="https://www.healthcare.siemens.com/ihe">https://www.healthcare.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="https://www.healthcare.siemens.com/hl7">https://www.healthcare.siemens.com/hl7</a>		
<b>DICOM</b>	<a href="https://www.healthcare.siemens.com/dicom">https://www.healthcare.siemens.com/dicom</a>		
<b>Link to General Information on IHE</b>			
<b>IHE</b>	<a href="https://www.ihe.net/IHE_Worldwide/">https://www.ihe.net/IHE_Worldwide/</a>		

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

---

Siemens Healthineers Headquarters  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
[www.siemens-healthineers.com](http://www.siemens-healthineers.com)

---

*Published by Siemens Healthineers · Doc-Ident. 11269744-AND-IHE · ©Siemens Healthineers CT, 2022*

---