

IHE Integration Statement

syngo CT VB10A

Product Name

NAEOTOM Alpha
SOMATOM go.All
SOMATOM go.Fit
SOMATOM go.Now
SOMATOM go.OpenPro
SOMATOM go.Sim
SOMATOM go.Top
SOMATOM go.Up
SOMATOM On.site
SOMATOM Pro.Pulse
SOMATOM X.ceed
SOMATOM X.cite

The Statement is also valid for regional product variants and configurations.



IHE Integration Statement		Date	10/2023
Vendor	Product Name	Version	
Siemens Healthineers	NAEOTOM Alpha, SOMATOM go.All, SOMATOM go.Fit, SOMATOM go.Now, SOMATOM go.OpenPro, SOMATOM go.Sim, SOMATOM go.Top, SOMATOM go.Up, SOMATOM On.site, SOMATOM Pro.Pulse, SOMATOM X.ceed, SOMATOM X.cite	VB10A	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF.b)	Acquisition Modality	Patient based worklist query Broad worklist query MPPS - Modality group case Siemens solution: Multi Study	
	Image Display	No options defined	
	Evidence Creator	none	
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined	
Radiation Exposure Monitoring (REM)	Acquisition Modality	none	
Access to Radiology Information (ARI)	Image Display	none	
Consistent Presentation of Images (CPI)	Print Composer	No options defined	
Consistent Time (CT)	Time Client	No options defined	
Portable Data for Imaging (PDI)	Portable Media Creator	No options defined	
	Image Display	No options defined	
	Print Composer	No options defined	
Internet address for vendor's IHE information: https://www.siemens-healthineers.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	https://www.siemens-healthineers.com/hl7		
DICOM	https://www.siemens-healthineers.com/dicom		
Link to General Information on IHE			
IHE	https://www.ihe.net/IHE_Worldwide/		

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
www.siemens-healthineers.com

Published by Siemens Healthineers · Doc-Ident. 11594610-AND-IHE · ©Siemens Healthineers CT, 2023
