Siemens Healthineers, P41 , ECM: XP_516141
11020808 ESK 023 05
Convert date: 2025-11-12T06:34:50 UTC
For signatures see info sheet (appended page)
Document is approved

IHE Integration Statement

MOBILETT Impact VC10/VC11



Published by Siemens Healthineers AG ·Doc-Ident. 11020808-ESK-023-05 *©Siemens Healthineers AG, 2025*

Siemens Healthineers, P41 , ECM: XP_516141
11020808 ESK 023 05
Convert date: 2025-11-12T06:34:50 UTC
For signatures see info sheet (appended page)
Document is approved

IHE Integration Statement		Date October 2025
Vendor Siemens Healthineers	**Product Name** MOBILETT Impact	**Version** VC10/VC11

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Acquisition Modality	Patient Based Worklist Query
		Assisted Acquisition Protocol Setting
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined
Consistent Presentation of Images (CPI)	Acquisition Modality	No options defined
	Print Composer	No options defined
Radiation Exposure Monitoring (REM)	Acquisition Modality	None

Internet address for vendor's IHE information: www.siemens-healthineers.com/ihe

Links to Standards Conformance Statements for the Implementation

HL7	www.siemens-healthineers.com/hl7
DICOM	www.siemens-healthineers.com/dicom

Links to general information on IHE

North America:	Europe:	Asia/Pacific:
ihe.net/IHE_Worldwide	ihe-europe.net	ihe.net/IHE_Worldwide

©Siemens Healthineers AG, 2025

Siemens Healthineers, P41 , ECM: XP_516141
11020808 ESK 023 05
Convert date: 2025-11-12T06:34:50 UTC
For signatures see info sheet (appended page)
Document is approved

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

Siemens Healthineers Headquarters
Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Published by Siemens Healthineers AG ·Doc-Ident. 11020808-ESK-023-05 *©Siemens Healthineers AG, 2025*

- generated automatically by SAP system **P41** -

Appendix to Document: **11020808 ESK 023 05 , ECM: XP_516141**
Sheet generated at : **2025-11-12T06:52:32 UTC**
Originator : **Siemens Healthineers, P41**
Signatures related to this document and performed in SAP:

Meaning	UTC date and time	surname, given name of signee
AUTHOR	**2025-11-04T15:20:02**	**Ke, Feng**
APPROVAL	**2025-11-12T06:51:59**	**Ni, Ming**
REVIEW	2025-11-07T01:44:51	Huang, Jing Hua
REVIEW	2025-11-12T06:45:53	Huang, Jing Hua